

# MWD Crossover Compact Electric Wheelchair

EIPW60




Customer				Ship To			
Account				Name			
Account #		Date		Address			
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #		City		State	
Buyer				ZIP	Phone		
Mark For				Email / Fax			

Please submit orders to: [magicmobility@sunmed.com](mailto:magicmobility@sunmed.com)

Seating Specifications			
HCPC	K0848   K0856   K0861   K0850   K0858   K0862 K0849   K0857   K0851   K0859	Typical WC Weight w/ Batteries	Starting at 320 lbs. Weights vary depending on options
Power Seating	50° of CG Tilt / Power Recline / 12" Seat Elevator	Seat Type	MPS and Rehab Style
Seat Size Range	Width: 13-23"   Depth: 13-24"	STF Height Range	17", 18", 19"
Drive Base Specifications			
Drive Wheel Position	Mid Wheel Drive	Drive Control	R-Net PM 120A controller
User Weight Capacity	350 lbs.	Drive Motors	4-pole 6mph
Base Size: Off-Road	Width 26"   Length 39"   Turning Radius 22.5"	Drive Wheels	14" Off-Road Drive Wheels   14" Urban Drive Wheels
Base Size: Urban	Width 24"   Length 39"   Turning Radius 22.5"	Tires / Casters	Pneumatic or Solid inserts
		Charger	Dual Mode
		Battery Size	Group 24
		Ground Clearance	3"
		Color Accents	10 colors

Client Measurements

- a) Top of head to seat \_\_\_\_\_ in
- b) Top of shoulder to seat \_\_\_\_\_ in
- c) Underarm to seat \_\_\_\_\_ in
- d) Elbow to seat \_\_\_\_\_ in
- e) Top of leg to seat \_\_\_\_\_ in
- f) Shoulder width \_\_\_\_\_ in
- g) Chest width \_\_\_\_\_ in
- i) Hip width \_\_\_\_\_ in
- j) Actual sitting depth \_\_\_\_\_ in
- m) Foot length \_\_\_\_\_ in
- n) Knee to foot \_\_\_\_\_ in
- o) Seat width \_\_\_\_\_ in



Notes

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
<b>Seating Options</b>			<b>\$</b>
<input type="checkbox"/> 600SF4		<b>REHAB</b> Seat and Back - Back Canes and adjustable flat seat pan	NC
<input type="checkbox"/> 600SF5		<b>MPS</b> Seat and Back - Captain style seating	NC
<input type="checkbox"/> 600SF6		<b>MPS</b> Back and <b>REHAB</b> Seat - Captain style backrest and adjustable flat seat pan	NC
<b>Base</b>			
<b>Magic 360 Base - Includes Transit and lights</b>		End Item Number: EIPW60	<b>\$</b>
<input type="checkbox"/> 600MW61	K0848	Magic Mobility - 360 No Power Seating Options	\$ 17,702.00
<input type="checkbox"/> 600MW62	K0856	Magic Mobility - 360 Single Power Seating Options	\$ 17,702.00
<input type="checkbox"/> 600MW63	K0861	Magic Mobility - 360 Multi Power Seating Options	\$ 17,702.00
<input type="checkbox"/> 600MW64	K0849	Magic Mobility - 360 No Power Seating Options - MPS (Captain)	\$ 17,702.00
<input type="checkbox"/> 600MW65	K0857	Magic Mobility - 360 Single Power Seating Options - MPS (Captain)	\$ 17,702.00
<input type="checkbox"/> 600MW66	K0862	Magic Mobility - 360 Multi Power Seating Options - MPS (Captain)	\$ 17,702.00
<b>Drive Wheels</b>			<b>\$</b>
<input type="checkbox"/> 600DW80		Off-Road Drive Wheels and Tires	NC
<input type="checkbox"/> 600DW81		Urban Drive Wheels and Tires - Pneumatic	NC
<input type="checkbox"/> 600DW82		Urban Drive Wheels and Tires - Solid	\$ 155.00
<b>Optional Drive Tire Options</b>			<b>\$</b>
<input type="checkbox"/> 6008A90		Tire Sealant for 2 Pneumatic Drive Tires	\$ 38.00
<b>Base Accessories</b>			<b>\$</b>
<input type="checkbox"/> 6008A91		Jack for Drive Tire	\$ 203.00
<b>Frame Color</b>			<b>\$</b>
<input type="checkbox"/> 600C185		Midnight Black	NC
<input type="checkbox"/> 600C186		Lucent Lime	NC
<input type="checkbox"/> 600C187		Sunset Red	NC
<input type="checkbox"/> 600C188		Horizon Blue	NC
<input type="checkbox"/> 600C189		Dusk Purple	NC
<input type="checkbox"/> 600C190		Vivid Pink	NC
<input type="checkbox"/> 600C191		Desert Orange	NC
<input type="checkbox"/> 600C192		Camouflage - Woods	NC
<input type="checkbox"/> 600C193		Camouflage - Desert	NC
<input type="checkbox"/> 600C171		Carbon Fiber	NC
<input type="checkbox"/> 600C194		American Flag	NC
<b>Seat To Floor Height</b>			<b>\$</b>
Seat to Floor Heights start at 17" +/- 0.25"			
<input type="checkbox"/> 600HT44		Lowest Seat To Floor Height	NC
<input type="checkbox"/> 600HT45		Seat To Floor Height +1"	NC
<input type="checkbox"/> 600HT46		Seat To Floor Height +2"	NC
<b>Batteries</b>			<b>\$</b>
<input type="checkbox"/> 600AG25	E2363	Group 24 Gel Battery (Pair) - Installed	\$ 1,116.00
<b>Battery Charger</b>			<b>\$</b>
<input type="checkbox"/> 600BJ28		10A Battery Charger	NC
<b>Power Retractable Docking Pin</b>			<b>\$</b>
<input type="checkbox"/> 6009A14		Retractable Docking Pin - Compatible with a Q-Straint QLK 150 docking system	\$ 3,071.00

**Seating: Seat Functions**

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
No Power Seat Options			\$
<input type="checkbox"/> 600ST87		No Power Seating System with Power Functionality	\$ 5,320.00
Power Tilt			\$

Utilized without Power Recline.

<input type="checkbox"/> 600PA64	E1002	50° Power Tilt in Space	\$ 8,292.00
----------------------------------	-------	-------------------------	-------------

Power Seat Elevate			\$
--------------------	--	--	----

Only available when Power Tilt or Power Tilt and Recline is selected.

<input type="checkbox"/> 600PI16	E2300	12" Adjustable Seat Elevator	\$ 4,750.00
----------------------------------	-------	------------------------------	-------------

Power Recline			\$
---------------	--	--	----

Utilized without Power Tilt. Requires No Power Seating with Power Functionality.

<input type="checkbox"/> 600PA65	E1003	150° Power Recline Only	\$ 6,750.00
----------------------------------	-------	-------------------------	-------------

Power Tilt and Recline			\$
------------------------	--	--	----

<input type="checkbox"/> 600PA66	E1006	50° Tilt & 150° Power Recline	\$ 12,750.00
----------------------------------	-------	-------------------------------	--------------

**Seating: Seat Sizes**

MPS (Captain) Seat Cushion (width x depth)			\$
--	--	--	----

MPS Backrest can be paired with the REHAB flat pan. If flat pan is desired, omit the MPS cushion and select the desired REHAB flat pan.

<input type="checkbox"/> 600SI24		16" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI25		18" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI26		18" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI27		20" x 18" Vinyl Seat	NC
<input type="checkbox"/> 600SI28		20" x 20" Vinyl Seat	NC
<input type="checkbox"/> 600SI29		22" x 20" Vinyl Seat	\$ 463.00

REHAB Flat Pan Seat Sizes - Adjustable in width and depth for most sizes			\$
--	--	--	----

**REHAB Seat Pan**

Standard seat pan available 13-23" wide. Deep seat pan available 13"-23" wide.

<input type="checkbox"/> 600PN13		Standard Adjustable Flat Seat Pan - 13" to 19" deep	NC
<input type="checkbox"/> 600PN14		Deep Adjustable Flat Seat Pan - 20" to 22" deep	\$ 463.00
<input type="checkbox"/> 600PN15		Deep Adjustable Flat Seat Pan - 23" to 24" deep	\$ 817.00

**REHAB Seat Width**

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

**REHAB Seat Depth**

<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"	<input type="checkbox"/> 21"	<input type="checkbox"/> 22"	<input type="checkbox"/> 23"	<input type="checkbox"/> 24"
------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------	------------------------------

REHAB Seating Cushion Options			\$
-------------------------------	--	--	----

Note cushion size will be determined by Rehab Seat size selected above

<input type="checkbox"/> 600CU59	E2607	JAY® Ion® (Sizes 16" - 21" Wide)	\$ 451.00
<input type="checkbox"/> 600CU60	E2608	JAY® Ion® (Size 22" Wide)	\$ 567.00
<input type="checkbox"/> 600CU61	E2607	JAY® Union® (Sizes 16" - 21" Wide)	\$ 510.00
<input type="checkbox"/> 600CU62	E2608	JAY® Union® (Size 22" Wide)	\$ 623.00
<input type="checkbox"/> 600CU64	E2622	JAY® Fusion® (Sizes 16" - 21" Wide)	\$ 504.00
<input type="checkbox"/> 600CU65	E2623	JAY® Fusion® (Size 22" Wide)	\$ 601.00

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
<b>Lateral Thigh Support Pads and Hardware (must pick pad and mount)</b>			<b>\$</b>
<b>Lateral Thigh Support Pads</b>			
<input type="checkbox"/> 600KF22	E0953	Lateral Thigh Support Pad Short - Right (5"x5")	\$ 144.00
<input type="checkbox"/> 600KF23	E0953	Lateral Thigh Support Pad Short - Left (5"x5")	\$ 144.00
<b>Lateral Thigh Support Mounts</b>			
<input type="checkbox"/> 600KF26	E1028	Lateral Thigh Removable Mount - Right	\$ 286.00
<input type="checkbox"/> 600KF27	E1028	Lateral Thigh Removable Mount - Left	\$ 286.00
<b>Seating: Seatback Options</b>			
<b>MPS Backrest Options - Comes standard with Manual Recline</b>			<b>\$</b>
<b>Backrest Width</b>			
<input type="checkbox"/> 600BN46		Standard Backrest Width - 16-20" Wide	NC
<input type="checkbox"/> 600BN47		Wide Backrest Width - 22-23" Wide	\$ 336.00
<b>Backrest Height</b>			
<input type="checkbox"/> 600B80		Standard Backrest Height - 22" Tall	NC
<input type="checkbox"/> 600B81		Tall Backrest Height - 26" Tall	\$ 558.00
<b>Backrest Material</b>			
<input type="checkbox"/> 600BU106		Cloth Backrest	NC
<input type="checkbox"/> 600BU107		Vinyl Backrest	\$ 304.00
<b>MPS Back Accessories</b>			<b>\$</b>
<input type="checkbox"/> 600PH33		Push Handle - Needed for O2 holder	\$ 159.00
<input type="checkbox"/> 6009A13		Accessory Mount - Needed for aftermarket headrest	\$ 121.00
<b>REHAB Seating Backrest Options</b>			<b>\$</b>
Backrest will have the width of the MPS seat or flat pan. Please choose your back height. Backrest will come with 2 point hardware with a 500 lb. weight capacity.			
<b>JAY J3 Back Contour</b>			
<input type="checkbox"/> 600JP01	E2313	J3® PA Back for 16"-20" Wide - 2.2" of lateral thoracic support	\$ 613.00
<input type="checkbox"/> 600JP02	E2614	J3® PA Back for 22" Wide - 2.2" of lateral thoracic support	\$ 910.00
<input type="checkbox"/> 600JP11	E2615	J3® PL Back for 16"-20" Wide - 3" of lateral thoracic support	\$ 715.00
<input type="checkbox"/> 600JP13	E2620	J3® PD Back for 16"-20" Wide (18" and 21" shell height only) - 6" of lateral thoracic support	\$ 858.00
<b>REHAB Seating Backrest Options</b>			<b>\$</b>
<b>JAY® J3® Back Shell Options and Height Ranges</b>			
<input type="checkbox"/> 600JR14		16.5" tall shell (Upper Thoracic Short)	NC
<input type="checkbox"/> 600JR15		18" tall shell (Upper Thoracic Med)	NC
<input type="checkbox"/> 600JR16		19.5" tall shell (Upper Thoracic Tall)	NC
<input type="checkbox"/> 600JR17		21" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR27		22.5" tall shell (Shoulder Height)	NC
<input type="checkbox"/> 600JR28		24" tall shell (Shoulder Height)	NC
<b>REHAB Flip Forward Backrest Option</b>			<b>\$</b>
<input type="checkbox"/> 600OD23		<b>REHAB</b> Flip Forward Backrest Option	\$ 296.00

**Seating: Seatback Options**

✓ OPTION	HCPCS	DESCRIPTION	RETAIL
JAY® Lateral Thoracic Swing-Away Supports (Requires Pad and Mount) - J3 back only			\$
<b>JAY® Lateral Thoracic Support Pads</b>			
<input type="checkbox"/> 600JW35	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Right	\$ 144.00
<input type="checkbox"/> 600JW38	E0956	Lateral Thoracic (L x H) 100 mm x 127 mm / 4" x 5" Pad - Left	\$ 144.00
<input type="checkbox"/> 600JW37	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Right	\$ 144.00
<input type="checkbox"/> 600JW40	E0956	Lateral Thoracic (L x H) 152 mm x 100 mm / 6" x 4" Pad - Left	\$ 144.00
<b>JAY® Lateral Thoracic Swing-Away Mount</b>			
<input type="checkbox"/> 600JW41	E1028	Lateral Thoracic Swing Away Bracket - Right	\$ 286.00
<input type="checkbox"/> 600JW42	E1028	Lateral Thoracic Swing Away Bracket - Left	\$ 286.00
JAY® Anterior Thoracic Support			\$
<input type="checkbox"/> 600ZR10	E0960	Bodypoint Monoflex Medium Chest Strap - Pad 3.5" x 18.5" (total length 49.5")	\$ 155.00
<input type="checkbox"/> 600ZR11	E0960	Bodypoint Monoflex Large Chest Strap - Pad 3.5" x 21.5" (total length 52.5")	\$ 155.00
<input type="checkbox"/> 600K10	E0960	JAY® Anterior Thoracic Support - Medium ATS8302	\$ 174.00
<input type="checkbox"/> 600K11	E0960	JAY® Anterior Thoracic Support - Large ATS8303	\$ 174.00
Pelvic Positioning Options - Select One			\$
<input type="checkbox"/> 600A848		Positioning Belt, Auto Style Buckle, 50 mm x 1600 mm / 2" x 63"	NC
<input type="checkbox"/> 600A849	E0978	Positioning Belt, Aircraft Buckle, Padded 50 mm x 1300 mm / 2" x 51"	\$ 86.00
<input type="checkbox"/> 600A850	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt	\$ 199.00
<input type="checkbox"/> 600A851	K0108	Bodypoint - Evoflex Pelvic Stabilizer Belt with Extension	\$ 220.00
Headrest Options - MPS Style			\$
<input type="checkbox"/> 600NU12	E0955	Cloth Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 411.00
<input type="checkbox"/> 600NU13	E0955	Vinyl Headrest - MPS seat back only (12" x 4.5" x 6.25")	\$ 411.00
Headrest Options - REHAB Style			\$
<b>Headrest Pad</b>			
<input type="checkbox"/> 600HE39	E0955	Whitmyer® Plush Headrest Pad 254 mm / 10"	\$ 234.00
<input type="checkbox"/> 600HE46	E0955	Whitmyer® Adjust-A-Plush Headrest Pad Medium (Adjustable from 7 to 16" Wide x 4" Tall)	\$ 360.00
<b>Headrest Mount</b>			
<input type="checkbox"/> 600MH33	E1028	Whitmyer® Linx Headrest Mount with Quick Release	\$ 188.00
<input type="checkbox"/> 600MH34	E1028	Whitmyer® Linx-2 Headrest Mount with Quick Release	\$ 229.00
Armrest Style Options			\$
<input type="checkbox"/> 600EG53		Flip-Up Removable Armrests (Post Style) - Attaches to base. Height adj 3" & angle adj 90°-155°	\$ 440.00
<input type="checkbox"/> 600EG54		Cantilever Armrests (Reclining Flip Back) - Attaches to back. Height adj 4" & angle adj 90°	\$ 884.00
<input type="checkbox"/> 600EG55		Cantilever 2 Point Quick Release Post - Attaches front of armrest to seat frame	\$ 1,173.00
Armrest Pad Options			\$
<input type="checkbox"/> 600AB80		4" x 16" Tapered Pads - Right	\$ 54.00
<input type="checkbox"/> 600AE80		4" x 16" Tapered Pads - Left	\$ 54.00
<input type="checkbox"/> 600AB81		10" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/> 600AE81		10" Waterfall Armrest Pad - Left	NC
<input type="checkbox"/> 600AB82		14" Waterfall Armrest Pad - Right	NC
<input type="checkbox"/> 600AE82		14" Waterfall Armrest Pad - Left	NC

Legrest Options			
<input checked="" type="checkbox"/> OPTION	HCPCS	DESCRIPTION	RETAIL
Center Mount Fixed Legrest			\$
Center Mount Fixed Legrest			
<input type="checkbox"/> 600H210	K0040	Center Mount Flip-Up Footplate	\$ 454.00
Center Mount Fixed Legrest - Calfpads			
<input type="checkbox"/> 600DI59	E0995	Fixed Center Mount Calf Pads	\$ 575.00
Center Mount Power Elevating Legrest			\$
Center Mount Power Elevating Legrest			
<input type="checkbox"/> 600H213	E1012	Center Mount Power Elevating - 2 Piece Footplate - Includes calf pads	\$ 4,175.00
Center Mount Power Elevating Legrest Extension			
<input type="checkbox"/> 600E42		Center Mount Power Elevating Extension - Standard - Footplate range 5" - 15"	NC
<input type="checkbox"/> 600E43		Center Mount Power Elevating Extension - Short - Footplate range 5" - 12.5"	NC
Center Mount Power Elevating Legrest Footplates			
<input type="checkbox"/> 600F117		Center Mount Power Elevating Footplate - Standard	NC
<input type="checkbox"/> 600F118		Center Mount Power Elevating Footplate - Short	NC
Swing-Away Legrest			\$
Swing-Away Legrest			
<input type="checkbox"/> 600H211	K0040	Side Mount Swing-Away Legrest - 2 piece flip-up footplate (8.5" x 6"). 13" height adjustable	\$ 795.00
S/A Non-Power Positioning Options			
<input type="checkbox"/> 600LS15	E0951	Heel and Side Support - Attaches to the footplate stem and footplate	\$ 657.00
<input type="checkbox"/> 600LS16	E0952	Foot Strap - 1.5" wide strap	\$ 89.00
Swing-Away Power Elevating Legrest			\$
Swing-Away Power Elevating Legrest			
<input type="checkbox"/> 600H214	E1010	Side Mount Power Elevating Legrest - Swing away legrest with flip-up 2" height adj footplate	\$ 4,555.00
S/A Power Legrest Positioning Options			
<input type="checkbox"/> 600DI58	K0042	2" Metal Heel Guards for 2 Piece Footplates - May not fit with all seat sizes	\$ 847.00
<input type="checkbox"/> 600HL15	E0951	Heel Loops - Attaches to the footplate. Length 11.5-13" and Width 1"	\$ 126.00
Accessories			\$
<input type="checkbox"/> 6008A93		Waterproof Wheelchair Cover	\$ 165.00
<input type="checkbox"/> 6008A94		Waterproof Joystick Cover	\$ 51.00
<input type="checkbox"/> 6008A95		Medical Necessity Bag - Velcro's to armrest (10.5" x 6" x 1.25")	\$ 38.00
<input type="checkbox"/> 600A221		Deluxe Backpack	\$ 108.00
<input type="checkbox"/> 600A615		Cup Holder - Self-leveling. Attaches to armrest opposite joystick	\$ 178.00
<input type="checkbox"/> 600A719	E2208	Articulating Oxygen Holder for <b>MPS</b> Seating System	\$ 1,171.00
<input type="checkbox"/> 600A612	E2207	Cane Holder for <b>MPS</b> Seating System - Includes push handles	\$ 753.00
<input type="checkbox"/> 6008A96		Luggage Rack for 360 - Attaches to base behind the chair (10.5" x 13")	\$ 374.00
<input type="checkbox"/> 6008A97		Charge Mate - 12 Volt/Car outlet and USB/cell phone charger with max draw of 6 amps	\$ 600.00
<input type="checkbox"/> 600A613		Cell Phone Holder - Adjustable to fit most cell phones	\$ 329.00
<input type="checkbox"/> 600A614		Fishing Rod Holder	\$ 329.00
Be Adaptive Products			
<input type="checkbox"/> 6008A98		Be Adaptive Limited Mobility Shooting Mount - Gun or crossbow	\$ 1,676.00
<input type="checkbox"/> 6008A99		Be Adaptive Trigger Mechanism	\$ 929.00
<input type="checkbox"/> 6008A100		Be Adaptive Limited Mobility Fishing Mount	\$ 1,676.00

Electronics			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
<b>Joystick</b>			\$
<input type="checkbox"/> 600J121		Standard LED Joystick (w/o Screen Display)	NC
<input type="checkbox"/> 600J116		Color Joystick (2.3" Color LCD Screen Display, CJSM1)	NC
<input type="checkbox"/> 600J118		High Resolution Color Joystick (CJSM2+BT) 3.5" Color LCD Display	\$ 1,206.00
<b>Joystick Mount</b>			\$
<input type="checkbox"/> 600JM34		Fixed on Armrest - Depth adjustable 8.25"	NC
<input type="checkbox"/> 600JM45	E1028	Swing-Away Joystick Mount - Includes latch to lock in the forward position	\$ 633.00
<input type="checkbox"/> 600JM46	E1028	Midline Joystick Mount - 2 position mount (Rotation 0° - 180°)	\$ 1,531.00
<b>Joystick Mount - Location</b>			\$
<input type="checkbox"/> 600HM16		Right Hand Armrest	NC
<input type="checkbox"/> 600HM17		Left Hand Armrest	NC
<b>Display Language Selection R-net</b>			\$
<input type="checkbox"/> 600L1		English	NC
<input type="checkbox"/> 600L3		French	NC
<input type="checkbox"/> 600L7		Spanish	NC
<b>Joystick Handle Type</b>			\$
Comes standard with PG Cone-Shaped Joystick Handle			
<input type="checkbox"/> 600JH13	E2323	Ball Handle	\$ 144.00
<input type="checkbox"/> 600JH26	E2323	Bodypoint Grooved Mushroom Joystick Handle	\$ 187.00
<input type="checkbox"/> 600JH24	E2323	Bodypoint – Rubber Dome Joystick Handle	\$ 187.00
<input type="checkbox"/> 600JH25	E2323	Bodypoint – U-Shaped Joystick Handle (4 inch)	\$ 187.00
<b>Expandable Power Module</b>			\$
Must select both items			
<input type="checkbox"/> 600EX05	E2377	Expandable 120A Controller - R-net	\$ 749.00
<input type="checkbox"/> 600EX06	E2313	Harness for 120A Expandable Controller	\$ 676.00
<b>Power Seat Function</b>			\$
<input type="checkbox"/> 600EO37	E2310	Thru Drive Control 1 Actuator - Must be selected if ordering 1 actuator function	\$ 2,196.00
<input type="checkbox"/> 600EO38	E2311	Thru Drive Control 2+ Actuators - Must be selected if ordering 2 or more actuator functions	\$ 3,975.00
<b>Power Seat Control Option</b>			\$
<input type="checkbox"/> 600QG11		Ctrl+5 Switch Box with Free Assignable Functions	\$ 1,430.00
<input type="checkbox"/> 600QG12		Ctrl+5 Switch Box with Free Assignable Functions (Multi-power only)	NC
<b>Ctrl+5 Switch Type</b>			\$
<input type="checkbox"/> 600EY22		Push Button	NC
<input type="checkbox"/> 600EY23	K0108	Toggle Switches	\$ 343.00
<b>Power Seat Control Option - Mounting location</b>			\$
<input type="checkbox"/> 600EY19		CTRL+5 MNT ARMREST RIGHT	NC
<input type="checkbox"/> 600EY20		CTRL+5 MNT ARMREST LEFT	NC
<b>Attendant Controls</b>			\$
<input type="checkbox"/> 6000A05	E2331	Dual Control Module for Attendant - Includes accessory mount and push handle	\$ 1,487.00
<b>Attendant Control Mounting Location</b>			\$
<input type="checkbox"/> 600DN5		Right Back Cane	NC
<input type="checkbox"/> 600DN6		Left Back Cane	NC



Electronics			
✓ OPTION	HCPCS	DESCRIPTION	RETAIL
R-net Omni-2 - 3.5" LCD Display w/ Built-In Infrared & Bluetooth			\$
<input type="checkbox"/> 600ON14		R-net Omni-2 Only Input Control Module	\$ 2,031.00
<input type="checkbox"/> 600ON15		R-net Omni2 Bundled Input Control Module Available w/ Switch-It Input Devices	NC
R-net Omni-2 -Mounting Location			\$
<input type="checkbox"/> 600DS23		Right Hand Armrest	NC
<input type="checkbox"/> 600DS24		Left Hand Armrest	NC
Head Control Packages			\$
<input type="checkbox"/> 600EU25	E2328	Switch-It DualPro® Proportional Head Array - Facial Pads	\$ 6,537.00
<input type="checkbox"/> 600EU26	E2328	Switch-It DualPro® Proportional Head Array - Spot Pads	\$ 6,537.00