



JAY[®] Zip[™] Cushion

E2607

January 2024

Mark For: _____

Date: _____

Dealer Acct #: _____

Dealer: _____

Dealer Contact: _____

Dealer Address: _____

Dealer City: _____ ST: _____ ZIP: _____

Dealer Phone: () _____ Fax: () _____

Confirmation Email: _____

Confirm Via: Fax Email

Submitting for: Quote Order

PO#: _____

ADDITIONAL SHIPPING INFORMATION

Ship To: _____

Attention: _____

Address: _____

Address: _____

Ship To City: _____ ST: _____ ZIP: _____

Ship To Phone: () _____ Fax: () _____

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.

JAY ZIP CUSHION - E2607

JAY Zip Cushion

<input type="checkbox"/>	JZXXYY	E2607	JAY Zip Cushion 8-16 wide	\$466.00
--------------------------	--------	-------	---------------------------	----------

Note: XX = Width and YY = Depth (Select appropriate assembly # on page 2)

Includes: Ionic+[™] fabric outer cover, inner cover, and foam base

2 Outer Covers included with each Cushion!




		JAY Zip Size Matrix									
		WIDTH									
DEPTH		8	9	10	11	12	13	14	15	16	
	8					N/A	N/A	N/A	N/A	N/A	
	9						N/A	N/A	N/A	N/A	
	10							N/A	N/A	N/A	
	11								N/A	N/A	
	12									N/A	
	13										
	14										
	15	N/A									
16	N/A	N/A									
17	N/A	N/A	N/A								
18	N/A	N/A	N/A	N/A							

	=	Standard Offering
JYW	=	Offered Via Jay Your Way
N/A	=	Not available

JAY Zip Measuring Guide

Measurements made from completed cushion with covers.



Bottom View

Width

X = Width

Note: Measurement made from bottom outside edge to outside edge of cushion.

Ex: JZ1214 X = 12"

Depth

Y = Depth

Note: Measurement made from bottom rear edge to front lower edge of cushion.

Ex: JZ1214 Y = 14"

Outer Cover

The JAY Zip outer cover is designed to promote comfort, convenience and breathability by the use of "Ionic+[™]" material and 3DX Spacer Fabric.



Inner Cover

The water resistant inner cover, combined with the Aquaguard[™] Zipper and anti wicking thread provides the ultimate protection against odor.



Note: To order via JAY Your Way add a "M" to the front of the part number and follow the prompts. For additional modifications see Page 2 or visit www.sunrisemedical.com to view even more!

Please visit www.sunparts.us for spare parts information.

JAY Zip Cushion

	Option #	HCPCS	Definition	Price		Option #	HCPCS	Definition	Price
<input type="checkbox"/>	JZ0808	E2607	JAY ZIP 8" X 8"	\$466.00	<input type="checkbox"/>	JZ1212	E2607	JAY ZIP 12" X 12"	\$466.00
<input type="checkbox"/>	JZ0809	E2607	JAY ZIP 8" X 9"	\$466.00	<input type="checkbox"/>	JZ1213	E2607	JAY ZIP 12" X 13"	\$466.00
<input type="checkbox"/>	JZ0810	E2607	JAY ZIP 8" X 10"	\$466.00	<input type="checkbox"/>	JZ1214	E2607	JAY ZIP 12" X 14"	\$466.00
<input type="checkbox"/>	JZ0811	E2607	JAY ZIP 8" X 11"	\$466.00	<input type="checkbox"/>	JZ1215	E2607	JAY ZIP 12" X 15"	\$466.00
<input type="checkbox"/>	JZ0812	E2607	JAY ZIP 8" X 12"	\$466.00	<input type="checkbox"/>	JZ1216	E2607	JAY ZIP 12" X 16"	\$466.00
<input type="checkbox"/>	JZ0813	E2607	JAY ZIP 8" X 13"	\$466.00	<input type="checkbox"/>	JZ1217	E2607	JAY ZIP 12" X 17"	\$466.00
<input type="checkbox"/>	JZ0814	E2607	JAY ZIP 8" X 14"	\$466.00	<input type="checkbox"/>	JZ1218	E2607	JAY ZIP 12" X 18"	\$466.00
<input type="checkbox"/>	JZ0908	E2607	JAY ZIP 9" X 8"	\$466.00	<input type="checkbox"/>	JZ1310	E2607	JAY ZIP 13" X 10"	\$466.00
<input type="checkbox"/>	JZ0909	E2607	JAY ZIP 9" X 9"	\$466.00	<input type="checkbox"/>	JZ1311	E2607	JAY ZIP 13" X 11"	\$466.00
<input type="checkbox"/>	JZ0910	E2607	JAY ZIP 9" X 10"	\$466.00	<input type="checkbox"/>	JZ1312	E2607	JAY ZIP 13" X 12"	\$466.00
<input type="checkbox"/>	JZ0911	E2607	JAY ZIP 9" X 11"	\$466.00	<input type="checkbox"/>	JZ1313	E2607	JAY ZIP 13" X 13"	\$466.00
<input type="checkbox"/>	JZ0912	E2607	JAY ZIP 9" X 12"	\$466.00	<input type="checkbox"/>	JZ1314	E2607	JAY ZIP 13" X 14"	\$466.00
<input type="checkbox"/>	JZ0913	E2607	JAY ZIP 9" X 13"	\$466.00	<input type="checkbox"/>	JZ1315	E2607	JAY ZIP 13" X 15"	\$466.00
<input type="checkbox"/>	JZ0914	E2607	JAY ZIP 9" X 14"	\$466.00	<input type="checkbox"/>	JZ1316	E2607	JAY ZIP 13" X 16"	\$466.00
<input type="checkbox"/>	JZ0915	E2607	JAY ZIP 9" X 15"	\$466.00	<input type="checkbox"/>	JZ1317	E2607	JAY ZIP 13" X 17"	\$466.00
<input type="checkbox"/>	JZ1008	E2607	JAY ZIP 10" X 8"	\$466.00	<input type="checkbox"/>	JZ1318	E2607	JAY ZIP 13" X 18"	\$466.00
<input type="checkbox"/>	JZ1009	E2607	JAY ZIP 10" X 9"	\$466.00	<input type="checkbox"/>	JZ1411	E2607	JAY ZIP 14" X 11"	\$466.00
<input type="checkbox"/>	JZ1010	E2607	JAY ZIP 10" X 10"	\$466.00	<input type="checkbox"/>	JZ1412	E2607	JAY ZIP 14" X 12"	\$466.00
<input type="checkbox"/>	JZ1011	E2607	JAY ZIP 10" X 11"	\$466.00	<input type="checkbox"/>	JZ1413	E2607	JAY ZIP 14" X 13"	\$466.00
<input type="checkbox"/>	JZ1012	E2607	JAY ZIP 10" X 12"	\$466.00	<input type="checkbox"/>	JZ1414	E2607	JAY ZIP 14" X 14"	\$466.00
<input type="checkbox"/>	JZ1013	E2607	JAY ZIP 10" X 13"	\$466.00	<input type="checkbox"/>	JZ1415	E2607	JAY ZIP 14" X 15"	\$466.00
<input type="checkbox"/>	JZ1014	E2607	JAY ZIP 10" X 14"	\$466.00	<input type="checkbox"/>	JZ1416	E2607	JAY ZIP 14" X 16"	\$466.00
<input type="checkbox"/>	JZ1015	E2607	JAY ZIP 10" X 15"	\$466.00	<input type="checkbox"/>	JZ1417	E2607	JAY ZIP 14" X 17"	\$466.00
<input type="checkbox"/>	JZ1016	E2607	JAY ZIP 10" X 16"	\$466.00	<input type="checkbox"/>	JZ1418	E2607	JAY ZIP 14" X 18"	\$466.00
<input type="checkbox"/>	JZ1108	E2607	JAY ZIP 11" X 8"	\$466.00	<input type="checkbox"/>	JZ1512	E2607	JAY ZIP 15" X 12"	\$466.00
<input type="checkbox"/>	JZ1109	E2607	JAY ZIP 11" X 9"	\$466.00	<input type="checkbox"/>	JZ1513	E2607	JAY ZIP 15" X 13"	\$466.00
<input type="checkbox"/>	JZ1110	E2607	JAY ZIP 11" X 10"	\$466.00	<input type="checkbox"/>	JZ1514	E2607	JAY ZIP 15" X 14"	\$466.00
<input type="checkbox"/>	JZ1111	E2607	JAY ZIP 11" X 11"	\$466.00	<input type="checkbox"/>	JZ1515	E2607	JAY ZIP 15" X 15"	\$466.00
<input type="checkbox"/>	JZ1112	E2607	JAY ZIP 11" X 12"	\$466.00	<input type="checkbox"/>	JZ1516	E2607	JAY ZIP 15" X 16"	\$466.00
<input type="checkbox"/>	JZ1113	E2607	JAY ZIP 11" X 13"	\$466.00	<input type="checkbox"/>	JZ1517	E2607	JAY ZIP 15" X 17"	\$466.00
<input type="checkbox"/>	JZ1114	E2607	JAY ZIP 11" X 14"	\$466.00	<input type="checkbox"/>	JZ1518	E2607	JAY ZIP 15" X 18"	\$466.00
<input type="checkbox"/>	JZ1115	E2607	JAY ZIP 11" X 15"	\$466.00	<input type="checkbox"/>	JZ1613	E2607	JAY ZIP 16" X 13"	\$466.00
<input type="checkbox"/>	JZ1116	E2607	JAY ZIP 11" X 16"	\$466.00	<input type="checkbox"/>	JZ1614	E2607	JAY ZIP 16" X 14"	\$466.00
<input type="checkbox"/>	JZ1117	E2607	JAY ZIP 11" X 17"	\$466.00	<input type="checkbox"/>	JZ1615	E2607	JAY ZIP 16" X 15"	\$466.00
<input type="checkbox"/>	JZ1209	E2607	JAY ZIP 12" X 9"	\$466.00	<input type="checkbox"/>	JZ1616	E2607	JAY ZIP 16" X 16"	\$466.00
<input type="checkbox"/>	JZ1210	E2607	JAY ZIP 12" X 10"	\$466.00	<input type="checkbox"/>	JZ1617	E2607	JAY ZIP 16" X 17"	\$466.00
<input type="checkbox"/>	JZ1211	E2607	JAY ZIP 12" X 11"	\$466.00	<input type="checkbox"/>	JZ1618	E2607	JAY ZIP 16" X 18"	\$466.00

JAY Your Way - Modifications

REDUCED HEIGHT			
<input type="checkbox"/>	RH	Remove 1/2" of foam from bottom	NC
NOTCH FOR BACK POST			
<input type="checkbox"/>	NP	Notch for Back Post	NC
Specify distance from back of cushion:			
		<input type="checkbox"/> 1" <input type="checkbox"/> 2" <input type="checkbox"/> 3"	
NOTCH FOR SEAT RAILS			
<input type="checkbox"/>	NR	Notch for Seat Rails	NC
COVER OPTIONS			
<input type="checkbox"/>	CVR	Cover Options	NC
		<input type="checkbox"/> Dartex <input type="checkbox"/> Incontinent <input type="checkbox"/> Microclimatic	

ANTERIOR OR POSTERIOR SLOPE			
<input type="checkbox"/>	WGP	Posterior Slope - Wedge	NC
Wedge glued underneath cushion (Cover sewn to match)			
Select desired height increase at front of cushion			
		<input type="checkbox"/> .5" <input type="checkbox"/> 1" <input type="checkbox"/> 1.5"	
<input type="checkbox"/>	WGA	Anterior Slope - Wedge	NC
Wedge glued underneath cushion (Cover sewn to match)			
Select desired height increase at rear of cushion			
		<input type="checkbox"/> .5" <input type="checkbox"/> 1" <input type="checkbox"/> 1.5"	

Don't see the option you want? NO PROBLEM! Write your request below and we can evaluate it and provide a quote.

Notes:



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA