



ORDER FORM

Squiggles + Stander

Prone, upright, and supine standing with up to 60° abduction.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Squiggles+ Stander

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Multi-Position Packages			\$
Package includes Squiggles+ Standard Model Multi-Position, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, and Small Sandals with Straps.			
<input type="checkbox"/> PSPSTB	Squiggles Plus Stander Package - Blue	E0641	\$6,766
<input type="checkbox"/> PSPSTP	Squiggles Plus Stander Package - Purple	E0641	\$6,766
<input type="checkbox"/> PSPSTR	Squiggles Plus Stander Package - Red	E0641	\$6,766
<input type="checkbox"/> PSPSTRV	Squiggles Plus Stander Package - Red Vinyl	E0641	\$6,766
Standard Model			\$

Selection Required.

Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplates and tray attachment.

<input type="checkbox"/>	184-600-P	Standard Support Shell - Prone	E0638	\$2,900
	184-764	Squiggles+ Pivot Base		\$1,552
<input type="checkbox"/>	184-600-S	Standard Support Shell - Supine	E0638	\$3,000
	184-764	Squiggles+ Pivot Base		\$1,552
<input type="checkbox"/>	184-600-M	Standard Support Shell - Multi-Position (Prone, Supine, and Upright)	E0641	\$3,173
	184-764	Squiggles+ Pivot Base		\$1,552
Positioning Pads			\$	

Selection Required. Includes chest and hip pad covers and knee support covers.

<input type="checkbox"/>	184-1800-02	Squiggles Plus Stander Positioning Pads - Blue		\$705
<input type="checkbox"/>	184-1800-04	Squiggles Plus Stander Positioning Pads - Purple		\$705
<input type="checkbox"/>	184-1800-09	Squiggles Plus Stander Positioning Pads - Red		\$705
<input type="checkbox"/>	184-1800-09VIN	Squiggles Plus Stander Positioning Pads - Red Vinyl		\$705

Accessories				\$
<input type="checkbox"/>	118-619	Tray - Gray		\$420
<input type="checkbox"/>	188-1854-05	Padded Tray Insert - Black		\$340
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray		\$305
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)		\$273
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)		\$298
<input type="checkbox"/>	117-787	Sandal Raiser 1" (Each) Requires Sandal Selection.	QTY	\$95

Headrest				\$
-----------------	--	--	--	-----------

Cushion selection required with hardware selection.

Headrest selection NOT available with Prone Configuration.

<input type="checkbox"/>	117-697	Headrest Support Hardware		\$236
		Flat Headrest		
<input type="checkbox"/>	184-1812-02	Flat Headrest Cushions - Blue		\$125
<input type="checkbox"/>	184-1812-04	Flat Headrest Cushions - Purple		\$125
<input type="checkbox"/>	184-1812-09	Flat Headrest Cushions - Red		\$125
<input type="checkbox"/>	184-1812-09VIN	Flat Headrest Cushions - Red Vinyl (Pair)		\$125
		Flat Headrest with Laterals		
<input type="checkbox"/>	184-1814-02	Flat Headrest and Lateral Cushions - Blue		\$218
<input type="checkbox"/>	184-1814-04	Flat Headrest and Lateral Cushions - Purple		\$218
<input type="checkbox"/>	184-1814-09	Flat Headrest and Lateral Cushions - Red		\$218
<input type="checkbox"/>	184-1814-09VIN	Flat Headrest and Lateral Cushions - Red Vinyl (Pair)		\$218

Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Distance Between Midline of Knees	5.5" - 8.4"	140 mm - 215 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.5" - 31.5"	480 mm - 800 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 mm L x 390 mm W x 390 mm H
Pivot Base Weight	14.3 lbs	6.5 kgs
Pivot Base Dimensions (Unfolded)	37.4" L x 23.2" W x 20.3" H	950 mm L x 590 mm W x 515 mm H
Pivot Base Dimensions (Folded)	26" L x 23.2" W x 20.3" H	660 mm L x 590 mm W x 515 mm H
Pivot Base Angle Range	Prone 90-70° / Supine 90-160°	Prone 90-70° / Supine 90-160°

Sunrise Medical (US) LLC
 2842 Business Park Avenue
 Fresno, CA 93727 USA
 Tel.: 1-800-333-4000
 Fax: 1-800-300-7502
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com



MK-129855 Rev. G ©05/2024



800.333.4000

www.SunriseMedical.com