



Empulse®



E|R90®

LET IT DRIVE YOU



red dot winner 2022
best of the best

PATENT-PENDING WHEEL UP

When manually propelling, raise the wheel for zero rolling resistance.

REMOVABLE BATTERY

Swap out the battery or remove it to make the R90 lighter to carry. A second or extended-range battery provides extra distance or even more convenience.

CENTER WHEEL POSITION

Positioned directly below the rear camber tube, the R90 minimally impacts the chair's footprint and center of gravity, providing unmatched maneuverability.

WIRELESS CONTROL

The wireless controller allows for a multitude of mounting positions and options to suit your individual needs.

SIMPLE, SECURE ATTACHMENT

The patent-pending V-shaped latch helps guide the R90 into its receiver and then lock it in place.

SUNRISE INTELLIGENCE APP

The app is available to track use, adjust max speed, and change drive profiles based on environment. In addition, the app customizes the drive performance to your weight and chair during the initial setup.



LEARN MORE



TECHNICAL DATA

HCPCS Code	E0896
Weight Capacity	275 lbs (125 kg)
Battery Range*	
Standard (2.5 Ah)	9.8 miles (15.7 km)
Extended (5 Ah)	19.5 miles (31.4 km)
Max Slope	6 Degrees /10.5%
Weight (without battery)	13.0 lbs (5.9 kg)
Standard Battery	1.5 lbs (0.7 kg)
Extended Battery	2.6 lbs (1.2 kg)
Frame Compatibility	Folding Frames Rigid Frames One-Arm Drive
Adjustable Max Speed	5.5 mph (9 kph)
Electronics	Bluetooth / Wireless
Battery	Lithium-Ion 36V 2.5Ah / 5Ah
Battery Charging Time	1 Hour (for 2.5Ah) 2 Hours (for 5Ah)
Motor	250 W / 22 Nm
Motor Brake	None

* Range is dependent on the ground surface, the weight of the rider, driving behavior, temperature and terrain.

WARNING

Using a power assist device can change your center-of-gravity or the way your chair rides. To reduce any risk of tip-over or fall, consult your healthcare provider about proper wheelchair configuration, including use of anti-tips and any other necessary adjustments.



This document can be viewed in PDF format at www.SunriseMedical.com.

SUNRISE MEDICAL CARES ABOUT YOUR SAFETY. Please read your owner's manual for safe operating instructions and maintenance guidelines.



Sunrise Medical (US), LLC, Fresno, CA, 93727
MK-130136 Rev. A ©02.2024



800.333.4000

www.SunriseMedical.com

