





ORDER FORM

Horizon Stander

EILK02

Robust and durable 3-in-1 standing frame for both young children and adult users.



Shipping Address						
Name						
Address						
City			S	tate		
ZIP		Phone				
Email / Fax						

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Horizon Stander			
Size 2 Weight Capaci	ty: 176 lbs Size 3 Weight Capacity: 220 lbs		
✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Frame			\$
Selection Required. I	ncludes Spine Assembly, Hip Pad, Chest Pad, Flexible Chest Laterals. (Knee supports and	Footplates separate)	
☐ 141FO5	Horizon Size 2 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	E0641	\$1,910
☐ 141FO6	Horizon Size 3 Extrusion with Chest Pad Hip pad Chest Laterals - Grey	E0641	\$1,990
Base			\$
Selection Required. I	MUST match Frame size. Includes worldwide adapters.		
☐ 141NB110	Size 2 Horizon Base/Chassis		\$3,528
☐ 141NB111	Size 3 Horizon Base/Chassis		\$4,403
Headrest			\$
Head Support Assem	bly selection required if Head Lateral Support/Extension Bracket is selected.		
☐ 141HE3	Size 2 and 3 Head Support Assembly - Black		\$287
☐ 141HE4	Size 2 and 3 Head Lateral Supports - Black (Pair)		\$248
☐ 141JG12	Headrest 4" Horizontal Extension Bracket		\$138
☐ 141Y029	141Y029 Size 1-3 Contoured Head Support		\$628
Hip Laterals			\$
Selection required if	Hip Belt is not selected.		
☐ 141KC7	Size 2 Horizon Hip Laterals - Black (Pair)	Size 2 Horizon Hip Laterals - Black (Pair)	
☐ 141KC8	Size 3 Horizon Hip Laterals - Black (Pair)		\$383
Hip Belts			\$
Selection required if	Hip Laterals are not selected. Multiple Selections can be made.		
☐ 141A856	Size 2 Supine Belt - Black	Size 2 Supine Belt - Black	
☐ 141A857	Size 3 Supine Belt - Black		\$386
Knee Supports			\$
	I Knee Assembly Required.		
Standard Knee Assen	nbly has 1"-3" adjustment or 1"-5" adjustment if flipped. Medium Knee Bracket has 5" - 7" a	djustment.	
Long Knee Bracket h			
☐ 141XX16	Size 2 and 3 Standard Knee Assembly - Black (Pair)		\$516
☐ 141JG14	Medium Knee Adjustment Bracket (Pair)		\$119
☐ 141JG15	Long Knee Adjustment Bracket (Pair)		\$122
Lower Leg Positionin			\$
Selection of Footplate			¢azz
☐ 1418A107	Small Sandals with Straps (Pair)		\$273 \$298
☐ 1418A108	Medium Sandals with Straps (Pair)		\$290
☐ 1418A109 ☐ 141F53		Large Sandals with Straps (Pair)	
☐ 141F65	Size 2 and 3 Filp up Footplate Assembly (Pair) Size 1-3 3D Footplate Assembly (Pair) - One size fits all	Size 2 and 3 Flip Up Footplate Assembly (Pair)	
Upper Extremity Supp	. , , , ,		\$1,332 \$
Must match Horizon			Ψ,
☐ 141A209	Size 2 and 3 Supine Tray		\$432
☐ 141A210	Size 2 and 3 Prone Tray		\$572
- ITIMEIU	SIZE 2 ditu 3 FTOTIE Hay		J J12

Technical Specifications - Imperial					
Frame Sizing Specification	Size 2	Size 3			
Age Range	9 18	14 Adult			
Horizontal Position L x W x H	50.4" × 25.2" × 31.7"	50.4" × 25.2" × 31.7"			
Weight of Product	110.6 lbs	130 lbs			
Maximum User Weight	176 lbs	220 lbs			
User Standing Height (Min / Max)	55.1" x 70.9"	61" / 70.9"			
Shoulder Height - Footplate to Top of Pad (Min / Max)	39.4" / 57"	39.4" / 57"			
Hip Pad Height - Footplate to Top of Pad (Min / Max)	29.5" / 37.4"	29.5" / 45.3"			
Hip Pad Width (Min / Max)	8.3" / 15.3"	12.2" / 18.5"			
Knee Height - Footplate to Middle of Pad (Min / Max)	15.7" / 21.6"	15.7" / 25.6"			
Width Between Chest Laterals (Min / Max)	9.4" / 14.9"	13.4" / 18.1"			
Footplate Angle Adjustable + / -	10°	10°			
Tray Sizes	20.5" x 23.2"	20.5" x 23.2"			
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°			
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°			
Frame Clearance - Bottom of Ski to Ground	7"	7"			

Technical Specifications - Metric					
Frame Sizing Specification	Size 2	Size 3			
Age Range	9 18	14 Adult			
Horizontal Position L x W x H	1280 mm x 640 mm x 805 mm	1280 mm x 640 mm x 805 mm			
Weight of Product	50 kgs	95 kgs			
Maximum User Weight	80 kgs	100 kgs			
User Standing Height (Min / Max)	1400 mm / 1800 mm	1550 mm / 1800 mm			
Shoulder Height - Footplate to Top of Pad (Min / Max)	1000 mm / 1450 mm	1000 mm / 1450 mm			
Hip Pad Height - Footplate to Top of Pad (Min / Max)	750 mm / 950 mm	750 mm / 1150 mm			
Hip Pad Width (Min / Max)	210 mm / 390 mm	310 mm / 470 mm			
Knee Height - Footplate to Middle of Pad (Min / Max)	400 mm / 550 mm	400 mm / 650 mm			
Width Between Chest Laterals [Min / Max]	240 mm / 380 mm	340 mm / 460 mm			
Footplate Angle Adjustable + / -	10°	10°			
Tray Sizes	520 mm x 590 mm	520 mm x 590 mm			
Tray Angle Adjustment - Prone Version	0° to 40°	0° to 40°			
Tray Angle Adjustment - Supine Version	-20° to +30°	-20° to +30°			
Frame Clearance - Bottom of Ski to Ground	180 mm	180 mm			

Sunrise Medical (US) LLC 2842 N. Business Park Avenue Fresno, CA 93727 USA



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com











