



ORDER FORM

# Squiggles Stander

The adjustable 3-in-1 stander designed for growing children.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

## Squiggles Stander

Weight Capacity: 48.4 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
<b>Packages</b>			\$
Package includes Squiggles Standard Model with Support Shell, Pivot Base, Positioning Pads, Flat Headrest Hardware, Cushion and Laterals, Tray, Small Sandals with Straps. Substitutions are not available for Packages.			
<input type="checkbox"/> PSSTB	Squiggles Stander Package - Blue	E0641	\$5,928
<input type="checkbox"/> PSSTBV	Squiggles Stander Package - Blue Vinyl	E0641	\$5,928
<input type="checkbox"/> PSSTG	Squiggles Stander Package - Green	E0641	\$5,928
<input type="checkbox"/> PSSTP	Squiggles Stander Package - Pink	E0641	\$5,928
<b>Standard Model</b>			\$

Selection Required.

Includes hip and chest support with sternum prompt; chest and hip lateral supports; size 2 knee cups, footplate and tray attachment.

<input type="checkbox"/> 118-600	Standard Model with Support Shell	E0641	\$2,218
<input type="checkbox"/> 118-764	Pivot Base - Pneumatic Angle Adjustment for 3-in-1 Positioning		\$1,620

<b>Positioning Pads</b>			\$
-------------------------	--	--	----

Selection Required. Includes chest and hip pad covers and knee support covers.


<input type="checkbox"/> 118-762	Squiggles Stander Positioning Pads - Blue		\$626
<input type="checkbox"/> 118-762-03VIN	Squiggles Stander Positioning Pads - Blue Vinyl		\$626
<input type="checkbox"/> 118-760	Squiggles Stander Positioning Pads - Green		\$626
<input type="checkbox"/> 118-763	Squiggles Stander Positioning Pads - Pink		\$626

<b>Accessories</b>			\$
--------------------	--	--	----

<input type="checkbox"/> 102-1640	QTY	Size 1 Knee Cup	\$128
<input type="checkbox"/> 102-3640	QTY	Size 3 Knee Cup	\$128
<input type="checkbox"/> 118-619		Tray - Gray	\$420
<input type="checkbox"/> 188-1854-05		Padded Tray Insert - Black	\$340
<input type="checkbox"/> 117-769		Grab Rail / Toy Bar for Tray	\$305
<input type="checkbox"/> 152-1600		Small Sandals w/ Straps (Pair)	\$273
<input type="checkbox"/> 152-2600		Medium Sandals w/ Straps (Pair)	\$298
<input type="checkbox"/> 120-816	QTY	Footplate Raiser 1" (Each)	\$92



Sandals



Grab Rail

## Squiggles Stander

✓ OPTION	DESCRIPTION	RETAIL
Headrest		\$

Cushion selection required with hardware.

Lateral Bracket Assembly Requires Flat Headrest Support Hardware and Cushion.

<input type="checkbox"/>	<b>120-859</b>	<b>Contoured Headrest Support Hardware</b>	\$385
<input type="checkbox"/>	137-800-03	Contoured Headrest Cushion - Blue	\$176
<input type="checkbox"/>	137-800-03VIN	Contoured Headrest Cushion - Blue Vinyl	\$176
<input type="checkbox"/>	137-800-01	Contoured Headrest Cushion - Green	\$176
<input type="checkbox"/>	137-800-04	Contoured Headrest Cushion - Pink	\$176
<input type="checkbox"/>	137-800-05	Contoured Headrest Cushion - Black	\$176
<input type="checkbox"/>	<b>117-697</b>	<b>Flat Headrest Support Hardware</b>	\$236
<input type="checkbox"/>	117-850-03	Flat Headrest Cushion - Blue	\$121
<input type="checkbox"/>	117-850-03VIN	Flat Headrest Cushion - Blue Vinyl	\$121
<input type="checkbox"/>	117-850-01	Flat Headrest Cushion - Green	\$121
<input type="checkbox"/>	117-850-04	Flat Headrest Cushion - Pink	\$121
<input type="checkbox"/>	117-850-05	Flat Headrest Cushion - Black	\$121
<input type="checkbox"/>	<b>118-667</b>	<b>Lateral Bracket Assembly Hardware (Pair)</b> Requires Flat Headrest Support Hardware and Cushion	\$217
<input type="checkbox"/>	137-807-03	Flat Headrest Lateral Cushions - Blue (Pair)	\$197
<input type="checkbox"/>	137-807-03VIN	Flat Headrest Lateral Cushions - Blue Vinyl (Pair)	\$197
<input type="checkbox"/>	137-807-01	Flat Headrest Lateral Cushions - Green (Pair)	\$197
<input type="checkbox"/>	137-807-04	Flat Headrest Lateral Cushions - Pink (Pair)	\$197
<input type="checkbox"/>	137-807-05	Flat Headrest Lateral Cushions - Black (Pair)	\$197



Headrest Support



Flat Headrest



Lateral Cushions



Contoured Headrest

## Technical Specifications

Frame Sizing Specification	Imperial	Metric
Age (approximate)	1 - 5 years	1 - 5 years
Maximum User Weight	48.4 lbs	22 kgs
User Standing Height	28.5" - 43.3"	730 mm - 1100 mm
Chest Width Between Laterals	6.3" - 9"	160 mm - 230 mm
Chest Support Height	21.7" - 32"	550 mm - 820 mm
Chest Support Angle Adjustment	+/- 30°	+/- 30°
Chest Support Depth Adjustment	2"	50 mm
Hip Width Between Hip Guides	6.3" - 9"	160 mm - 230 mm
Hip Support Height from Footplate to Mid Buttocks	11.8" - 20.8"	300 mm - 530 mm
Footplate Angle	10° Plantarflexion / 10° Dorsiflexion	10° Plantarflexion / 10° Dorsiflexion
Small Sandal Size	6.3" L x 3" W	175 mm L x 75 mm W
Medium Sandal Size	8.3" L x 3.2" W	210 mm L x 80 mm W
Tray Height: Prone	21.9" - 32.5"	555 mm - 825 mm
Tray Height: Supine	18.9" - 32.5"	480 mm - 825 mm
Tray Angle Adjustment	Prone 30° / Supine 40°	Prone 30° / Supine 40°
Stander Support Frame Weight	15.4 lbs	7 kgs
Stander Support Frame Dimensions	26.8" L x 15.4" W x 15.4" H	680 x 390 x 390 mm
Pivot Base Weight	12.1 lbs	5.5 kgs
Pivot Base Dimensions (Unfolded)	32.3" L x 21.7" W x 20.3" H	820 x 550 x 515 mm
Pivot Base Dimensions (Folded)	27.6" L x 21.7" W x 11.8" H	700 x 550 x 300 mm
Knee Cup Opening Width	Size 1: 2.5"      Size 3: 3.5"	Size 1: 63 mm      Size 3: 89 mm

### Sunrise Medical (US) LLC

2842 Business Park Avenue  
 Fresno, CA 93727 USA  
 Tel.: 1-800-333-4000  
 Fax: 1-800-300-7502  
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at [www.SunriseMedical.com](http://www.SunriseMedical.com)



MK-129857 Rev. F ©01/2024



800.333.4000

[www.SunriseMedical.com](http://www.SunriseMedical.com)