

Q200R


QUICKIE®




Q200R

Ultra-Transportable, Everyday Performance.

Q200R

The perfect balance of size and performance.



BALANCE

With the Q200 R, our vision was to create a true indoor/outdoor power chair, one with super-tight indoor maneuverability and ZERO compromise on outdoor stability and performance.

RELIABLE

The chassis is formed from a single sheet of metal, and its steel components receive automotive undercoat plating. The industry's best motors undergo additional end-of-line testing before leaving our factory.

PERFORMANCE

The Q200 R provides great traction when driving up ramps because the chair weight shifts towards the rear drive wheels. When traveling outside at higher speeds, the Q200 R tracks straight so you are in complete control.



PERFECT FOR VEHICLE LIFTS: Starting at just 216 lbs., the Q200 R works great with hoist and boom style vehicle lifts. Worried about space? With a fully fold-down back and removable leg supports/arm supports, it'll comfortably fit in most SUVs or vans.



Captain's Seat

The Captain Seat provides a relaxed sitting position and features a semi-recline fold-down back, adjustable seat depth, and convertible seat pan that allows you to use your preferred cushion. Available with a 22° tilt and power elevating, swing away leg supports.



Rehab Seat

The Rehab Seat is an adjustable seat frame with a sleek look. Its minimal design alleviates weight capacity concerns when using aftermarket equipment. With two seat sizes 16" and 19" the arm supports can be adjusted to accommodate 15" to 21" widths.

Q200 R SEAT FRAME COMPARISON

	CAPTAIN SEAT	REHAB SEAT
CryoFluid™ (Available only on JAY® Fusion)	•	•
Manual Recline Option	•	
Power 22° Center of Gravity Tilt	•	•
Power Leg Support	•	•
Seat Width	16", 18", 20"	15" to 21"
Seat Depth	16" to 22"	16" to 20"
Back Height	25"	18" to 30"


Q200 R

See order form as specifications may vary based upon configuration

HCPC codes	No Power Seating Options	Single Power Seating Options
Group 3 Standard Duty	K0848	K0856
Max. user weight	300 lbs.	
Seat-to-floor height	17" to 19"	
Total width	22.8" *	
Total length (without footplates)	32"	
Min. turning radius	18.3"	
Max. speed	6 mph	
Battery size	Group 22	
Overall weight	Starting at 216 lbs.	
Drive wheel position	Rear	
Max. ground clearance	3"	
Wheel options	12" or 14" drive wheels, 9" casters	
Max. range	Up to 15.5 miles	
Transit securement brackets	ISO 7176-19:2008 Transit Approved	
Warranty	Lifetime Limited	

* Exact measurements may vary depending on options.

COLOR OPTIONS



Important Note—
Print processes cannot accurately represent painted colors. Actual colors may vary.



Order Form



Web Page



Sunrise Medical (US) LLC, Fresno, CA 93727
MK-129983 Rev.C ©05.2022

800.333.4000 www.SunriseMedical.com

