

Pulse™







QUICKIE PUISE

Available in 6 Models:

Quickie Pulse 5

BC (K0848)

5 MPH, 2-Pole motor, No Power Option base with ASAP II Rehab Seat

CC (K0849)

5 MPH, 2-Pole motor, No Power Option base with Captain's Seat

Quickie Pulse 6

BC (K0848)

6 MPH, 4-Pole motor, No Power Option base with ASAP II Rehab Seat

CC (K0849)

6 MPH, 4-Pole motor, No Power Option base with Captain's Seat

SC (K0856)

6 MPH, 4-Pole motor, Single Power Option Tilt (SPOT) or Power Recline with ASAP II Rehab Seat or 3.7 Recline Seat

MPC (K0861)

6 MPH, 4-Pole motor, Power Tilt and Recline with 3.7 Recline Seat



- Mid wheel drive
- Reliable and compact for indoor accessibility
- Flexible suspension system for optimal outdoor performance
- Black or gray tire color options
- Available with 15 shroud color options
- Reliable 2-pole and 4-pole motors simulated a mean time of 20 years to failure during High Cycle Fatigue test

Power Seating Options

- NEW 3.7 Power Recline with Power Shear Reduction
- 50° of Power Tilt to optimize positioning
- 9" Power Seat Elevate for increased access to cabinets and eye-to-eye interaction with peers

Simple to use VR2 non-expandable electronics with 4-way toggle switch for operation of power seating without the need for expandable controls

Power Elevating Legrests with Articulation

Expandable and Non-Expandable Controls

Electronics Options

between environments







ASL Micro Joystick





Specialty input devices available on-chair and will ship complete with Pulse Series wheelchairs

R-net expandable electronics with assignable buttons for control of common functions and 8 programmable profiles for transitions

- 16 Switch-It® proportional and non-proportional (switch) input control options available
- 8 ASL® proportional and non-proportional (switch) input control options available



Switch-It Head Array



Multiple Seat Frame Options

Captain's Seat (CC)

- Fully adjustable seat size, legrest placement, and armrest width
- Convertible seat pan for flexibility in cushion selection
- Semi-recline and fold-down backrest
- Optional 02, crutch, and walker holders

ASAP II Rehab Seat (BC & SC)

- In-field width, depth, armrest, legrest, and back height adjustability
- Manual reclining and fold-down backrest
- Available with back upholstery, SPOT Curved backrest, JAY® & ART SPO customizable backrests, and the JAY J3 Back
- Extensive seating and positioning capabilities
- Optional 02, crutch, and walker holders

3.7 Recline Rehab Seat (SC & MPC)

- Option of Recline Curved back or Integrated JAY J3 Back for seamless functionality with power recline and shear reduction system
- Low profile electronics enclosure shell for a clean look and compatibility with accessories
- Seat-to-floor heights as low as 16.3"
- Cantilever Reclining Armrests stay parallel to ground when reclined

NEW Options for Consumer Convenience

- RAM® Phone and Tablet Mounts allow mobile phones and/or tablets to be connected and disconnected for easy viewing
- Active Controls Power Buddy provides an auxiliary USB power source for charging devices and optional LED lights for improved visibility
- RAM® Mounts Self-Leveling Cup Holder uses gimbal system to keep drinks level during travel or when stopped on inclines
- . Active Controls ReCharge™ with Socket eliminates the challenges of connecting and disconnecting the battery charger by providing an easy-to-reach adapter



Phone Mount



Tablet Mount



Power Buddy



ReCharge™



Quickie Pulse Specifications

SPECIFICATIONS	PULSE 5 BC	PULSE 6 BC	PULSE 5 CC	PULSE 6 CC	PULSE 6 SC	PULSE 6 MPC
Drive Wheel Position	Mid	Mid	Mid	Mid	Mid	Mid
HCPCS Code	K0848	K0848	K0849	K0849	K0856	K0861
Speed	5.0 MPH	6.0 MPH	5.0 MPH	6.0 MPH	6.0 MPH	6.0 MPH
Overall Weight	235 lbs.	235 lbs.	250 lbs.	250 lbs.	260 lbs.	345 lbs.
Overall Width	24"	24"	24"	24"	24"	24"
Overall Length	34"	34"	34"	34"	34"	34"
User Weight Capacity	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.	300 lbs.
Seat Widths	12" - 22"	12" to 20"	16", 18", 20"	16", 18", 20"	12" to 22"	14" to 22"
Seat Depths	12" - 22"	12" to 22"	15" to 22"	15" to 22"	12" to 22"	16" to 20"
Seat To Floor Height	16" - 20"	16" to 20""	16" to 18" (from pan)	16" to 18" (from pan)	16.3" to 21.6"	16.3" to 18.3"
Seating Frame Options	ASAP II Rehab Seat	ASAP II Rehab Seat	Captain's Seat	Captain's Seat	ASAP II Rehab Seat or 3.7 Recline	3.7 Recline
Transit Approved	ISO	ISO	WC-19*	WC-19*	Tilt: ISO/Recline: WC-19*	WC-19*
Warranty	Limited/Lifetime	Limited/Lifetime	Limited/Lifetime	Limited/Lifetime	Limited/Lifetime	Limited/Lifetime

^{*} Conforms with ANSI/RESNA WC-19 standards. Please refer to product Owner's Manual for additional information. Owner's Manuals are available for download at www.sunrisemedical.com





