



ORDER FORM

Mygo Seat

Versatile therapy chair that grows with its user without sacrificing proximal pelvic positioning and support.

Customer			
Account			
Account #		Date	
<input type="checkbox"/> Quote	<input type="checkbox"/> Order	PO #	
Buyer			
Mark For			

Shipping Address			
Name			
Address			
City		State	
ZIP		Phone	
Email / Fax			

The HCPCS CODES herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative HCPCS CODES, including misc. codes to ensure access for their enrollees. The use of HCPCS CODES does not ensure coverage or payment.

Mygo Seat

Size 1 Weight Capacity: 110 lbs | Size 2 Weight Capacity: 132 lbs

✓ OPTION	DESCRIPTION	HCPCS	RETAIL
Seat Shell			\$
Selection Required.			
<input type="checkbox"/> 117-610S	Size 1 Mygo Seat Shell		
117-610-1	Seat Base with Depth Adjustment	E2294	\$911
117-610-2	Seat Back - Fixed	E2293	\$911
117-610-3	Seat Back Reclining Mechanism	E1014	\$479
117-610-4	Seat Base Hip Guides / Support	E0956	\$133
117-610-5	Seat Base Medial Thigh Guides	E0957	\$166
117-610-6	Seat Base Leg Adductors	K0108	\$133
<input type="checkbox"/> 137-600S	Size 2 Mygo Seat Shell		
137-600-1	Seat Base with Depth Adjustment	E2294	\$1,011
137-600-2	Seat Back - Fixed	E2293	\$1,011
137-600-3	Seat Back Reclining Mechanism	E1014	\$479
137-600-4	Seat Base Hip Guides / Support	E0956	\$133
137-600-5	Seat Base Medial Thigh Guides	E0957	\$166
137-600-6	Seat Base Leg Adductors	K0108	\$133
Seat Cover			\$
Selection Required. Must match Seat Shell size.			
<input type="checkbox"/> 117-1851-08	Size 1 Mygo Seat Cover - Blue		N/C
<input type="checkbox"/> 117-1851-06	Size 1 Mygo Seat Cover - Gray		N/C
<input type="checkbox"/> 117-1851-07	Size 1 Mygo Seat Cover - Orange		N/C
<input type="checkbox"/> 117-1851-09	Size 1 Mygo Seat Cover - Pink		N/C
<input type="checkbox"/> 117-2851-08	Size 2 Mygo Seat Cover - Blue		N/C
<input type="checkbox"/> 117-2851-06	Size 2 Mygo Seat Cover - Gray		N/C
<input type="checkbox"/> 117-2851-07	Size 2 Mygo Seat Cover - Orange		N/C
<input type="checkbox"/> 117-2851-09	Size 2 Mygo Seat Cover - Pink		N/C
Base			\$
<input type="checkbox"/> 102-2753	Size 2 Hi Low Chassis (Foot Pedal Operated. Includes Push Handle)	E1399	\$2,251
<input type="checkbox"/> 194-700	Mobility Base Interface Plate (Allows mounting Mygo Seat onto Sunrise Wheelchairs)	K0108	\$525
Pelvic Fixation			\$
Selection of either Pelvic Harness or Pelvic Cradle required.			
Hip Laterals require Pelvic Cradle			
<input type="checkbox"/> 117-801	Pelvic Harness & Brackets Small - Black	E0978	\$164
<input type="checkbox"/> 117-802	Pelvic Harness & Brackets Medium - Black	E0978	\$164
<input type="checkbox"/> 120-789	Pelvic Harness Spacer Pads - Black (Pair)	K0108	\$92
<input type="checkbox"/> 117-477	Pelvic Cradle Size 3 (old size 0)	K0108	\$454
<input type="checkbox"/> 117-478	Pelvic Cradle Size 4 (old size 1)	K0108	\$504
<input type="checkbox"/> 137-632	Hip Laterals (Pair)	E0956	\$338
Dynamic Backrest Kit			\$
<input type="checkbox"/> 137-661	Dynamic Backrest 100N Cylinder	K0108	\$466
Armrest			\$
<input type="checkbox"/> 117-1777	Size 1 Mygo Adjustable Armrests (Pair)	E0973	\$332
<input type="checkbox"/> 117-2777	Size 2 Mygo Adjustable Armrests (Pair)	E0973	\$365

Mygo Seat

✓ OPTION	DESCRIPTION	HCPCS	RETAIL	
Lower Leg Positioning			\$	
Footplate selection required with Sandals.				
Sandal Selection Required with Sandal Raisers.				
<input type="checkbox"/>	117-1776	Size 1 Mygo Angled Footplates (Pair)	K0040	\$332
<input type="checkbox"/>	117-2776	Size 2 Mygo Angled Footplates (Pair)	K0040	\$365
<input type="checkbox"/>	152-1600	Small Sandals w/ Straps (Pair)	E0951, E0952	\$273
<input type="checkbox"/>	152-2600	Medium Sandals w/ Straps (Pair)	E0951, E0952	\$298
<input type="checkbox"/>	152-3600	Large Sandals w/ Straps (Pair)	E0951, E0952	\$321
<input type="checkbox"/>	117-787	QTY Sandal Raiser 1" (Each)	K0108	\$95
<input type="checkbox"/>	117-765-05VIN	Mygo Knee Pads (vinyl)	K0108	\$227
Harness			\$	
<input type="checkbox"/>	155-N001-L600	Mygo/SQSE Trunk Harness - Small	E0960	\$406
<input type="checkbox"/>	155-N001-L601	Mygo/SQSE Trunk Harness - Medium	E0960	\$406
<input type="checkbox"/>	155-N001-L602	Mygo/SQSE Trunk Harness - Large	E0960	\$406
Lateral Support			\$	
Selection of Chest Support Requires Laterals. Not available with Complex Laterals.				
<input type="checkbox"/>	117-837	Flip-Away Laterals with Cover (Pair)	E0956, E1028	\$579
<input type="checkbox"/>	137-634	Flip-Away Laterals Extended Height with Cover (Pair) <i>Required w/ Size 2 Seat</i>	E0956, E1028	\$711
<input type="checkbox"/>	117-803	Rigid Laterals with Cover (Pair)	E0956	\$331
<input type="checkbox"/>	137-633	Rigid Laterals Extended Height with Cover (Pair)	E0956	\$465
<input type="checkbox"/>	137-637	Flip-Away Complex Lateral Horizontal with Cover Left	E0956, E1028	\$371
<input type="checkbox"/>	137-638	Flip-Away Complex Lateral Horizontal with Cover Right	E0956, E1028	\$371
<input type="checkbox"/>	137-808-05	Chest Support - Black	E0960	\$187
Headrest			\$	
Selection of Headrest Support Hardware (Contoured/Flat) Requires Cushion				
Headrest Laterals only available with Flat Headrest				
<input type="checkbox"/>	120-859	Contoured Headrest Support Hardware	E1028	\$385
<input type="checkbox"/>	137-800-05	Contoured Headrest Cushion - Black	E0955	\$176
<input type="checkbox"/>	117-697	Flat Headrest Support Hardware	E1028	\$236
<input type="checkbox"/>	117-850-05	Flat Headrest Cushion - Black	E0955	\$121
<input type="checkbox"/>	120-893-05	Flat Headrest Laterals with Covers - Black (Pair)	E0956	\$338
Tray			\$	
Armrest selection required with Tray				
<input type="checkbox"/>	117-646	Mygo Size 1 and 2 Black Tray	E0950	\$398
<input type="checkbox"/>	117-652-BLACK	Black Padded Tray (cover with tray)	K0108	\$621
<input type="checkbox"/>	117-769	Grab Rail / Toy Bar for Tray	K0108	\$305
Accessories			\$	
Shoulder Support is standard on Mygo Size 2.				
<input type="checkbox"/>	137-621	Shoulder Support Hardware for Mygo Size 1 Cushion selection required.		\$281
<input type="checkbox"/>	137-801-08	Shoulder Support Cushion - Blue		\$160
<input type="checkbox"/>	137-801-06	Shoulder Support Cushion - Gray		\$160
<input type="checkbox"/>	137-801-07	Shoulder Support Cushion - Orange		\$160
<input type="checkbox"/>	137-801-09	Shoulder Support Cushion - Pink		\$160

Technical Specifications - Universal

Specification	Size 1	Size 2
Age (approximate)	3 - 10 years	8 - 14 years
Backrest Angle: Prone	10°	10°
Backrest Angle: Recline	25°	25°
Footplate: Abduction	8°	8°
Footplate: Adduction	12°	20°
Footplate: Plantar/Dorsiflexion	10° / 10°	10° / 10°

Size 2 Hi Lo Chassis (Excluding Push Bar)

Tilt in Space: Prone	10°	10°
Tilt in Space: Recline	25°	25°

Technical Specifications - Imperial

Specification	Size 1	Size 2
User Height	41.3" - 59"	50" - 66"
User Weight	39.6 - 110 lbs	39.6 - 132 lbs
Seat Depth	10.6" - 16.5"	13.8" - 18.5"
Seat Width	8" - 12.8"	8.7" - 13.6"
Backrest Height	14.2' - 18.5"	18.1" - 22.4"
Seat to Sandal	8.5" - 13.8"	12.4" - 18.5"
Chest Width (Distance between laterals)	6.7" - 10.6"	6.7" - 10.6"
Armrest Height	6.3" - 8.3"	8.3" - 10.2"
Knee Width	3.5" - 4.3"	4.7" - 5.5"
Tray Size	21.6" x 18.9"	21.6" x 18.9"
Seat Unit Weight	22 lbs	31.9 lbs

Size 2 Hi Lo Chassis (Excluding Push Bar)

Top of Seat to Floor	13.4" - 25.8"	13.4" - 25.8"
Base Weight	26.5 lbs	26.5 lbs

Technical Specifications - Metric

Specification	Size 1	Size 2
User Height	1050 mm - 1500 mm	1270 mm - 1680 mm
User Weight	18 - 50 kg	18 - 60 kg
Seat Depth	270 mm - 420 mm	350 mm - 470 mm
Seat Width	200 mm - 325 mm	220 mm - 345 mm
Backrest Height	360 mm - 470 mm	460 mm - 570 mm
Seat to Sandal	215 mm - 350 mm	315 mm - 470 mm
Chest Width (Distance between laterals)	170 mm - 270 mm	170 mm - 270 mm
Armrest Height	160 mm - 210 mm	210 mm - 260 mm
Knee Width	90 mm - 110 mm	120 mm - 140 mm
Tray Size	550 mm x 480 mm	550 mm x 480 mm
Seat Unit Weight	10 kgs	14.5 kgs

Size 2 Hi Lo Chassis (Excluding Push Bar)

Top of Seat to Floor	340 mm - 655 mm	340 mm - 655 mm
Base Weight	12 kgs	12 kgs

Sunrise Medical (US) LLC

2842 Business Park Avenue
 Fresno, CA 93727 USA
 Tel: 1-800-333-4000
 Fax: 1-800-300-7502
 E-Mail: webmaster@sunmed.com



If you are visually impaired, this document can be viewed in PDF format at www.SunriseMedical.com



MK-129849 Rev. H ©01/2024



800.333.4000

www.SunriseMedical.com