



# Encompass Backrest

E2617

January 2024

**Mark For:** \_\_\_\_\_  
 Date: \_\_\_\_\_  
 Dealer Acct #: \_\_\_\_\_  
 Dealer: \_\_\_\_\_  
 Dealer Contact: \_\_\_\_\_  
 Dealer Address: \_\_\_\_\_  
 Dealer City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Dealer Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_  
 Confirmation Email: \_\_\_\_\_  
 Confirm Via:  Fax  Email

**Submitting for:**  Quote  Order  
 PO#: \_\_\_\_\_

**ADDITIONAL SHIPPING INFORMATION**

Ship To: \_\_\_\_\_  
 Attention: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Address: \_\_\_\_\_  
 Ship To City: \_\_\_\_\_ ST: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 Ship To Phone: ( ) \_\_\_\_\_ Fax: ( ) \_\_\_\_\_

The **HCPCS CODES** herein are based on PDAC verification or interpretation of Medicare definitions and guidelines. Non-Medicare payers may accept alternative **HCPCS CODES**, including misc. codes to ensure access for their enrollees. The use of **HCPCS CODES** does not ensure coverage or payment.


## STEP 1 - Select Seat Back

Encompass Backrest				
✓	Option #	HCPCS	Definition	Price
<input type="checkbox"/>	<b>SBK90XX</b>	<b>E2617</b>	<b>Solid Back Package</b>	<b>\$583</b>
<i>Select Size</i>				
<input type="checkbox"/>	ENC1000		18" Tall or Less	\$288
<input type="checkbox"/>	ENC1010		Greater than 18" Tall	\$288
	ENC1500		Removable Encompass Mounting Hardware	\$201
	HR0001		Universal Headrest Adapter Plate	\$94

Back sizes available in 1" increments only. Standard foam is 1" High Resiliency. Cover is black Air Exchange/contact surface; Dartex/non-contact surface. No additional upholstery options are available.

Width Dimensions				
<i>Wheelchair width determines correct selection of back, which is sized to fit between uprights.</i>				
<input type="checkbox"/>	14"	<input type="checkbox"/>	15"	<input type="checkbox"/>
<input type="checkbox"/>	16"	<input type="checkbox"/>	17"	<input type="checkbox"/>
<input type="checkbox"/>	18"	<input type="checkbox"/>	19"	<input type="checkbox"/>
<input type="checkbox"/>	20"	<input type="checkbox"/>	21"	<input type="checkbox"/>
<input type="checkbox"/>	22"	<input type="checkbox"/>	23"	<input type="checkbox"/>
<input type="checkbox"/>	24"	<input type="checkbox"/>	25"	<input type="checkbox"/>

Height Dimensions				
<input type="checkbox"/>	12"	<input type="checkbox"/>	13"	<input type="checkbox"/>
<input type="checkbox"/>	14"	<input type="checkbox"/>	15"	<input type="checkbox"/>
<input type="checkbox"/>	16"	<input type="checkbox"/>	17"	<input type="checkbox"/>
<input type="checkbox"/>	18"	<input type="checkbox"/>	19"	<input type="checkbox"/>
<input type="checkbox"/>	20"	<input type="checkbox"/>	21"	<input type="checkbox"/>
<input type="checkbox"/>	22"	<input type="checkbox"/>	23"	<input type="checkbox"/>
<input type="checkbox"/>	24"	<input type="checkbox"/>	25"	<input type="checkbox"/>
<input type="checkbox"/>	26"	<input type="checkbox"/>	27"	<input type="checkbox"/>
<input type="checkbox"/>	28"			



Shown without Foam or Cover

Encompass Sizing Guide						
Wheelchair Width	Lateral Size(s) (Standard)	Overall Back Width (w/o Laterals)	Overall Back Width with Posterior Laterals		Patient Thoracic Width Using Swing-Away Laterals	
			Minimum	Maximum	Minimum	Maximum
14	Small	7.2	8.5	14	7.5	13.5
15	Small	8.2	9.5	15	8.5	14.5
16	Small	9.2	10.5	16	9.5	15.5
17	Med-Small	10.2	11.5	17	10.5	16.5
18	Med-Small	11.2	12.5	18	11.5	17.5
19	Med-Small	12.2	13.5	19	12.5	18.5
20	Med-Large	13.2	14.5	20	13.5	19.5
21	Med-Large	14.2	15.5	21	14.5	20.5
22	Med-Large	15.2	16.5	22	15.5	21.5
23	Large	16.2	17.5	23	16.5	22.5
24	Large	17.2	18.5	24	17.5	23.5
25	Large	18.2	19.5	25	18.5	24.5

## STEP 2 - Select Foam Modification

Foam Modifications to Back section are in addition to the 1/2" closed cell base layer. Recommended total foam thickness is 2" including base layer.

T-Foam			
<input type="checkbox"/>	ENCFOAMT1	1" T-Foam	\$94
<input type="checkbox"/>	ENCFOAMT1.5	1.5" T-Foam	\$106
<input type="checkbox"/>	ENCFOAMT2	2" T-Foam	\$124
<input type="checkbox"/> X-Soft (Yellow) <input type="checkbox"/> Soft (Pink) <input type="checkbox"/> Medium (Blue)			

Gel			
<input type="checkbox"/>	ENCFOAMG.5	1/2" Gel	\$201
<input type="checkbox"/>	ENCFOAMG1	1" Gel	\$201
<input type="checkbox"/>	ENCFOAMGC	1/2" Cubed Gel	\$201

Sunmate Foam			
<input type="checkbox"/>	ENCFOAMSM1	1" Sunmate Foam	\$77
<input type="checkbox"/>	ENCFOAMSM1.5	1.5" Sunmate Foam	\$94
<input type="checkbox"/>	ENCFOAMSM2	2" Sunmate Foam	\$112
<input type="checkbox"/> Soft <input type="checkbox"/> Med. Soft <input type="checkbox"/> Medium <input type="checkbox"/> Firm			

Additional Foam			
<input type="checkbox"/>	ENCFOAMX.5	Add HR Foam per 1/2"	\$28
			Inches X \$25 =
<input type="checkbox"/>	ENC5000	F.I.P. Kit	\$499

<b>SUBTOTAL OF SECTION 2: FOAM MODIFICATION</b>
\$ _____

## STEP 3 - Select Laterals

Foam Modifications Laterals are in addition to the 1/2" closed cell base layer. Recommended total foam thickness is 1" including base layer.

### STEP 3A - Posterior Laterals

A. Posterior (POS) Laterals - LEFT ONLY				
<input type="checkbox"/>	ENC2010L-RT	E0956	LEFT SM 9.5"H X 3.25"D	\$146
<input type="checkbox"/>	ENC2020L-RT	E0956	LEFT MED-SM 10.5"H X 3.5"D	\$146
<input type="checkbox"/>	ENC2030L-RT	E0956	LEFT MED-LG 11.25"H X 3.5"D	\$146
<input type="checkbox"/>	ENC2040L-RT	E0956	LEFT LG 11.75"H X 3.5"D	\$146
<input type="checkbox"/>	ENC2099L-RT	E0956	LEFT Personalized POS Lateral	QUOTE
		Height	Depth	

A. Posterior (POS) Laterals - RIGHT ONLY				
<input type="checkbox"/>	ENC2010R-RT	E0956	RIGHT SM 9.5"H X 3.25"D	\$146
<input type="checkbox"/>	ENC2020R-RT	E0956	RIGHT MED-SM 10.5"H X 3.5"D	\$146
<input type="checkbox"/>	ENC2030R-RT	E0956	RIGHT MED-LG 11.25"H X 3.5"D	\$146
<input type="checkbox"/>	ENC2040R-RT	E0956	RIGHT LG 11.75"H X 3.5"D	\$146
<input type="checkbox"/>	ENC2099R-RT	E0956	RIGHT Personalized POS Lateral	QUOTE
		Height	Depth	

Notes:

Notes:

<b>SUBTOTAL OF SECTION 3A: POS LATERALS</b>
\$ _____

### STEP 3B - Swing Away Laterals

B. Swing-Away Anterior (S/A) Laterals - LEFT ONLY				
<i>Must select POS Lateral, same size (SM - LG) recommended</i>				
<input type="checkbox"/>	ENC2510L-RT	E1028	LEFT SM 6.25"H X 3.75"D	\$287
<input type="checkbox"/>	ENC2520L-RT	E1028	LEFT MED-SM 6.75"H X 4.25"D	\$287
<input type="checkbox"/>	ENC2530L-RT	E1028	LEFT MED-LG 7.25"H X 4.75"D	\$287
<input type="checkbox"/>	ENC2540L-RT	E1028	LEFT LG 7.75"H X 5.25"D	\$287
<input type="checkbox"/>	ENC2599L-RT	E1028	LEFT Personalized POS Lateral	QUOTE
		Height	Depth	

B. Swing-Away Anterior (S/A) Laterals - RIGHT ONLY				
<i>Must select POS Lateral, same size (SM - LG) recommended</i>				
<input type="checkbox"/>	ENC2510R-RT	E1028	RIGHT SM 6.25"H X 3.75"D	\$287
<input type="checkbox"/>	ENC2520R-RT	E1028	RIGHT MED-SM 6.75"H X 4.25"D	\$287
<input type="checkbox"/>	ENC2530R-RT	E1028	RIGHT MED-LG 7.25"H X 4.75"D	\$287
<input type="checkbox"/>	ENC2540R-RT	E1028	RIGHT LG 7.75"H X 5.25"D	\$287
<input type="checkbox"/>	ENC2599R-RT	E1028	RIGHT Personalized POS Lateral	QUOTE
		Height	Depth	

Notes:

Notes:

<b>SUBTOTAL OF SECTION 3B: S/A LATERALS</b>
\$ _____

## STEP 4 - Calculate Component Prices

### 4A. Base Package Price (Section 1+2)

*HCPCS - E2617*

Back Package (SBK90XX)	\$583
<b>ENTER Foam Mod. (Subtotal - Step 2)</b>	\$
<b>TOTAL</b>	\$

### 4B. Laterals Price (Section 3A+3B)

<b>ENTER POS Lateral (Subtotal - 3A)</b>	\$
<b>ENTER S/A Lateral (Subtotal - 3B)</b>	\$
<b>TOTAL</b>	\$

### GRAND TOTAL

4A. Back Package Price	\$
4B. Laterals Price	\$
<b>GRAND TOTAL</b>	\$

**NOTES**



Sunrise Medical (US) LLC  
2842 Business Park Ave. · Fresno, CA 93727 · USA