



QUICKIE P222 SE BD/SD

Solid Seat (No Power or Single Power Seating Options)

EIPW2 - Effective January 2019

Account: _____	Ship To: _____
Date: _____	Name: _____
PO#: _____	Address: _____
Buyer: _____	City, State, Zip: _____
Marked For: _____	Phone: _____

*The HCPCS CODES provided are based on code verification through PDAC or our interpretation of the code definitions.
The use of HCPCS CODES does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.*

✓	Option #	HCPCS	Definition	Price	✓	Option #	HCPCS	Definition	Price
---	----------	-------	------------	-------	---	----------	-------	------------	-------

1. QUICKIE P222 SE GROUP 4 CATEGORY

Quickie P222 SE				
<input type="checkbox"/>	116PC19	K0868 K0014	Quickie P222 SE BD (Solid Seat - No Power)	\$10,275.00
<input type="checkbox"/>	116PC1	K0877 K0014	Quickie P222 SE SD (solid seat/single power seat option)	\$10,275.00
<input type="checkbox"/>	116PC3	K0877 K0014	Quickie P222 SE (Base Only)	\$9,575.00

Performance Package				
<input type="checkbox"/>	116PF8	-	8.5 (8.5 mph / 350 LB Capacity) (Max. Wt. cap. For CG Tilt 250lbs)	N/C

Power CG Tilt (53° Tilt)				
<i>See electronic section for control options</i>				
<input type="checkbox"/>	116ST2	E1002	P222 SE Power CG Tilt Max. Wt. capacity for CG Tilt 250lbs. Not available with base only purchase.	\$5,150.00

Tilt Weight Capacity				
<input type="checkbox"/>	116LB3	-	Up to 175 lb Weight Capacity	N/C
<input type="checkbox"/>	116LB4	-	176 to 250 lb Weight Capacity	N/C

LED Light Package (Rnet only)				
<input type="checkbox"/>	116LI04	-	LED Light Package (NA Base only)	\$1,030.00

2. FRAME COLOR OPTIONS

Standard Colors				
<input type="checkbox"/>	116C140	-	Aztec Gold	NC
<input type="checkbox"/>	116C2	-	Black (Not avail as base w/Zebra)	NC
<input type="checkbox"/>	116C141	-	Black Cherry	NC
<input type="checkbox"/>	116C10	-	Black Opal	NC
<input type="checkbox"/>	116C38	-	Blue Opal	NC
<input type="checkbox"/>	116C36	-	Candy Blue	NC
<input type="checkbox"/>	116C14	-	Candy Purple	NC
<input type="checkbox"/>	116C35	-	Candy Red	NC
<input type="checkbox"/>	116C142	-	Electric Blue	NC
<input type="checkbox"/>	116C143	-	Evergreen	NC
<input type="checkbox"/>	116C126	-	Glow	NC
<input type="checkbox"/>	116C144	-	Green Apple	NC
<input type="checkbox"/>	116C147	-	Hot Sparkle Pink	NC
<input type="checkbox"/>	116C148	-	Mauve Pink	NC

<input type="checkbox"/>	116C103	-	Rootbeer	NC
<input type="checkbox"/>	116C114	-	Sparkle Silver	NC
<input type="checkbox"/>	116C145	-	Sunrise Orange	NC
<input type="checkbox"/>	116C113	-	Titanium Color Paint	NC
<input type="checkbox"/>	116C8	-	Yellow	NC

Standard Matte Colors				
<input type="checkbox"/>	116C116	-	Matte Black	NC
<input type="checkbox"/>	116C151	-	Matte Black Cherry	NC
<input type="checkbox"/>	116C152	-	Matte Purple	NC
<input type="checkbox"/>	116C150	-	Matte Electric Blue	NC
<input type="checkbox"/>	116C149	-	Matte Evergreen	NC

Custom Colors				
<input type="checkbox"/>	116C43	-	Mossy Oak Camouflage	\$305.00
<input type="checkbox"/>	116C42	-	Stars and Stripes	\$305.00
<input type="checkbox"/>	116C112	-	Desert Camo	\$305.00
<input type="checkbox"/>	116C108	-	Zebra (Choose base color above: White Std, n/a black)	\$305.00

3. BATTERIES and CHARGERS

Batteries (Gel Cells)				
<input type="checkbox"/>	116AG03	E2363	GP24 MK Gel Flat Top Batteries Pair w/strap	\$975.00
<input type="checkbox"/>	116AG22	E2363	GP24 Brightway AGM Flat Top Batteries Pair w/strap	\$735.00
<input type="checkbox"/>	116AG05	-	Omit batteries Includes harness & straps for standard post batteries	N/C

Battery Charger				
<input type="checkbox"/>	116BJ19	-	8A Dual Mode (Compact) CTE	N/C
<input type="checkbox"/>	116BJ3	-	Omit Charger	-\$40.00

4. SEATING

Seat Width				
<input type="checkbox"/>	116S5	-	14" (not avail. w/ 19" & 20" seat depths)	N/C
<input type="checkbox"/>	116S6	-	15" (not avail. w/ 19" & 20" seat depths)	N/C
<input type="checkbox"/>	116S7	-	16" (not avail. w/19" & 20" seat depths and dual post arms)	N/C
<input type="checkbox"/>	116S8	-	17"	N/C
<input type="checkbox"/>	116S9	-	18"	N/C
<input type="checkbox"/>	116S10	-	19"	N/C
<input type="checkbox"/>	116S11	-	20"	N/C
<input type="checkbox"/>	116S27	-	22" (not avail. w/CG Power Tilt; Base width is 26.5" w/ 12"x4" tires)	N/C
<input type="checkbox"/>	116S28	-	24" (not avail. w/CG Power Tilt; Base width is 26.5" w/ 12"x4" tires)	N/C

Seat Height				
<input type="checkbox"/>	116HT10	-	18" (not avail. w/long ext. tubes; add 1.75" for CG Tilt)	N/C
<input type="checkbox"/>	116HT12	-	19" (add 1.75" for CG Tilt)	N/C
<input type="checkbox"/>	116HT14	-	20" (add 1.75" for CG Tilt)	N/C

Seat Base Type				
<input type="checkbox"/>	116SL2	-	Seat Sling Parapak (n/a w/CG Power Tilt)	N/C
<input type="checkbox"/>	116SL5	-	Seat Pan Aluminum	N/C
<input type="checkbox"/>	116SL6	-	Omit Seat Sling	-\$15.00

Seat Depth				
<input type="checkbox"/>	116SD7	-	14"	N/C
<input type="checkbox"/>	116SD8	-	15"	N/C
<input type="checkbox"/>	116SD9	-	16"	N/C
<input type="checkbox"/>	116SD10	-	17"	N/C
<input type="checkbox"/>	116SD11	-	18"	N/C
<input type="checkbox"/>	116SD12	-	19"	N/C
<input type="checkbox"/>	116SD73	-	20"	N/C

5. SEAT CUSHIONS

Standard Cushion				
<input type="checkbox"/>	116CU1	K0669	2" Cushion.	\$120.00

JAY Cushions (Optional)			
<i>To order any JAY cushions on-chair, please see the JAY price list.</i>			
<i>Cushion name or part number must be filled out w/ desired width & depth</i>			
<input type="checkbox"/>	116JC20	JAY Bundled Cushion	Price List
		Cushion Name or Part #	
		Width: _____	Depth: _____

6. BACKS

Backrest Type				
<input type="checkbox"/>	116BT74	-	Folding Angle Adj. Backrest	\$115.00
<input type="checkbox"/>	116BT3	-	Non-Folding Angle Adj. Backrest	N/C

Back Height (Upholstery)				
<input type="checkbox"/>	116B4	-	Low 15.5"	N/C
<input type="checkbox"/>	116B8	-	Medium 17.5"	N/C
<input type="checkbox"/>	116B13	-	Tall 19.5"	N/C

Back Upholstery				
<input type="checkbox"/>	116BU27	-	Parapak Upholstery (select back height)	N/C
<input type="checkbox"/>	116BU26	-	Jay Solid Backrest (select back height; Decrease seat depth by 1")	N/C
<input type="checkbox"/>	116BU14	-	Omit Back Upholstery	-\$15.00

Jay Solid Backrest Back Height				
<input type="checkbox"/>	116B36	-	16"	N/C
<input type="checkbox"/>	116B37	-	17"	N/C
<input type="checkbox"/>	116B38	-	18"	N/C
<input type="checkbox"/>	116B39	-	19"	N/C
<input type="checkbox"/>	116B40	-	20"	N/C
<input type="checkbox"/>	116B41	-	21"	N/C
<input type="checkbox"/>	116B42	-	22"	N/C
<input type="checkbox"/>	116B43	-	23"	N/C
<input type="checkbox"/>	116B44	-	24"	N/C

Back Post Type				
<input type="checkbox"/>	116BP1	-	Straight Back Post	N/C
<input type="checkbox"/>	116BP41	-	8 Bend Back Post - Not avail. w/ Jay Solid backrest	\$105.00

Jay 3 Option				
<input type="checkbox"/>	116JB01	-	Jay J3 (select OMIT Back Uph above)	See J3 OF

7. POSTURAL SUPPORTS

Headrest (Incl. on Specified Specialty Controls)

JAY Solid Backrest Only

<input type="checkbox"/>	Headrest (8" x 4.5")		\$460.00
	116HE65	E0955 Curved Headrest 8" x 4.5"	\$230.00
	116MH1	E1028 Removable Mount	\$230.00
<input type="checkbox"/>	Headrest (9" x 4.5")		\$460.00
	116HE66	E0955 Curved Headrest 9" x 4.5"	\$230.00
	116MH1	E1028 Removable Mount	\$230.00

Pad Style

JAY Solid Backrest Only

<input type="checkbox"/>	Lateral Pads Curved		\$460.00
	116KE1	E0956 Lateral Pad Curved Pair	\$235.00
	116LM1	E1028 Swing Away Mount Pair	\$225.00
<input type="checkbox"/>	Lateral Pads Straight		\$450.00
	116KE2	E0956 Lateral Pad Straight Pair	\$225.00
	116LM1	E1028 Swing Away Mount Pair	\$225.00

Pad Size

<input type="checkbox"/>	116KD1	-	Lateral Support 3" x 6" (pr)	N/C
<input type="checkbox"/>	116KD2	-	Lateral Support 4" x 5" (pr)	N/C
<input type="checkbox"/>	116KD3	-	Lateral Support 4" x 6" (pr)	N/C
<input type="checkbox"/>	116KD4	-	Omit Lateral Support	-\$25.00

8. ARMRESTS

Armrest (Pair)

Single Post Height Adj. Armrests are N/C on Power CG Tilt

<input type="checkbox"/>	116AA131	-	Dual Post Flip Back Short (14-16in)	N/C
<input type="checkbox"/>	116AA132	-	Dual Post Flip Back Long (16-20in)	N/C
<input type="checkbox"/>	116AA51	E0973	Dual Post Hgt Adj - Short (14-16 in)	\$285.00
<input type="checkbox"/>	116AA52	E0973	Dual Post Hgt Adj - Long (16-20 in)	\$285.00
<input type="checkbox"/>	116AA8	E0973	Single Post Height Adjustable	\$270.00
<input type="checkbox"/>	116AA14	-	Omit Armrests and mounts	N/C

Armpads (Pair)

<input type="checkbox"/>	116AM7	-	Desk Length Waterfall (10") (n/a. w/ sgl. post Armrest w/Jay Solid Backrest)	N/C
<input type="checkbox"/>	116AM8	-	Full Length Waterfall (14") (n/a. w/ sgl. post Armrest w/Jay Solid Backrest)	N/C
<input type="checkbox"/>	116AM1	-	Desk Length Standard (10")	N/C
<input type="checkbox"/>	116AM2	-	Full Length Standard (14")	N/C

9. FOOTRESTS

Hanger (Pair)

<input type="checkbox"/>	116H12	-	70° Fixed Straight (only avail. w/ solid 1-piece footplate.; n/a. w/ CG Power Tilt)	N/C
<input type="checkbox"/>	116H2	-	70° Swing Away Hanger (Not. Avail w/ solid 1-piece footplate.)	N/C
<input type="checkbox"/>	116H13	-	80° Fixed Tapered (only avail. w/ solid 1-piece footplate.; n/a. w/ CG Power Tilt)	N/C
<input type="checkbox"/>	116H37	-	90° Swing Away (Only available with kids angle adjustable footplates)	N/C
<input type="checkbox"/>	116H16	E0990	Elevating Hanger (Req. aluminum footplate or angle adj.)	\$280.00
<input type="checkbox"/>	116H32	-	Omit Hangers	-\$15.00

Extension Tube

Not available with 90° hangers.

<input type="checkbox"/>	116E1	-	Short (13" - 15")	N/C
<input type="checkbox"/>	116E2	-	Medium (15" - 18") Not available with multi-position footplate.	N/C
<input type="checkbox"/>	116E5	-	Long (18" - 20") Not available with 70° or 80° hangers.	N/C

Footrest (Pair)

<input type="checkbox"/>	Composite w/ Heel Loops		N/C	
	Not Avail. w/ fixed hangers			
	116F1	-	Composite Footplate	N/C
	116HL1	E0951	Heel Loop	\$30.00
<input type="checkbox"/>	Adult Angle Adjust. w/ heel loop		\$220.00	
	Not Avail. w/ fixed hangers			
	116F282	K0040	Adult Angle Adj. Footplate	\$170.00
	116HL10	E0951	Heel Loops	\$50.00
<input type="checkbox"/>	Kids Angle Adjust. w/ heel loop		\$205.00	
	Only avail. w/ 90° hangers			
	116F283	K0040	Kids Angle Adj. Footplate	\$170.00
	116LS13	K0038	Leg Straps	\$35.00
<input type="checkbox"/>	116F287	-	Aluminum (Only avail. with ELRs)	\$35.00
<input type="checkbox"/>	116F60	K0040	Solid One Piece not avail. w/Power CG Tilt	\$35.00
<input type="checkbox"/>	Two-piece Locking w/Heel Loop		\$220.00	
	116F58	K0040	Two-piece Locking Footplate	\$170.00
	116HL10	E0951	Heel Loop	\$50.00
<input type="checkbox"/>	Multi-Positional w/Heel Loop		\$225.00	
	Not Avail. w/ fixed hangers			
	116F47	K0040	Multi-Positional Footplate	\$175.00
	116HL10	E0951	Heel Loop	\$50.00
<input type="checkbox"/>	116F6	-	Omit Foot Plate w/Heel Loop	N/C

Footrest Accessories

<input type="checkbox"/>	116LS3	-	Extra Leg Strap, Black	\$35.00
--------------------------	--------	---	------------------------	---------

10. ELECTRONICS

Drive Control Options

<input type="checkbox"/>	116ER04	R-Net Non-Expandable Drive Control (120A)	NC
<input type="checkbox"/>	R-Net Expandable Drive Control Upgrade (120A)		\$980.00
	116EX03	E2377 Expandable 120A Controller - R-net	\$540.00
	116EX04	E2313 Harnessing for 120A Expandable Controller	\$440.00

R-net Joystick (for expandable electronics)

<input type="checkbox"/>	116J70	-	Std LED Joystick for Expandable Ctrl's	N/C
<input type="checkbox"/>	116J108	-	Joystick w/Stereo Jacks & Color Display	\$210.00
<input type="checkbox"/>	116J88	-	Joystick W/Toggle, Spd Pot, Stereo Jacks & Color Display	\$410.00
<input type="checkbox"/>	116J95	-	Omit JS for use with OMNI	N/C

Joystick Mount (Fixed or Swing Away)

<input type="checkbox"/>	116JM10	-	Fixed Mount	N/C
<input type="checkbox"/>	116JM11	E1028	Swing Away Mount	\$245.00

Joystick Mount (Location)

Height Adjustable.

<input type="checkbox"/>	116HM1	-	Right Side Mount	N/C
<input type="checkbox"/>	116HM2	-	Left Side Mount	N/C

Joystick Handle (Optional)

<input type="checkbox"/>	116JH1	E2323	Ball Handle	\$80.00
<input type="checkbox"/>	116JH2	E2323	Foam Ball Handle	\$80.00
<input type="checkbox"/>	116JH4	E2323	T-Handle	\$80.00
<input type="checkbox"/>	116JH6	E2323	Goalpost	\$80.00

R-net OMNI and Output Options

(For expandable controls, choose all that apply)

<input type="checkbox"/>	116ON1	-	Rnet OMNI Input Control Module and display with built in Infrared and armrest mount. Choose mount below. (required for Specialty Input Devices.)	\$1,005.00
--------------------------	--------	---	--	------------

OMNI Mount

(If used with Joystick must be mounted opposite side)

<input type="checkbox"/>	116DS03	E1028	Swing Away Left Side OMNI Mount	\$245.00
<input type="checkbox"/>	116DS04	E1028	Swing Away Right Side OMNI Mount	\$245.00

Switch inputs for Joysticks or OMNI

Drive or Mode / Profile Control

<input type="checkbox"/>	116SW22	-	Stereo to 2 Mono Splitter	\$65.00
<input type="checkbox"/>	116SW6	-	Ribbon Switch	\$100.00
<input type="checkbox"/>	116SW1	-	Buddy Button Switch	\$80.00
<input type="checkbox"/>	116SW2	-	Disc Switch	\$165.00
<input type="checkbox"/>	116SW9	-	Wobble Switch	\$515.00
<input type="checkbox"/>	116SW4	-	Micro Light Switch	\$80.00
<input type="checkbox"/>	116SW24	-	Egg Switch-Black	\$100.00
<input type="checkbox"/>	116SW21	-	Dual toggle switch only	\$175.00
<input type="checkbox"/>	116SW23	-	4-way toggle switch w/5th button DB9 Connector.	\$175.00

Switch Mount Bracket

<input type="checkbox"/>	116EM01	-	Single Switch Mount Bracket	\$200.00
<input type="checkbox"/>	116EM02	-	Dual toggle switch mount bracket	\$200.00
<input type="checkbox"/>	116EM03	-	4-way toggle switch mount bracket.	\$200.00

Specialty Driver Control Options

Requires expandable controls and OMNI. All packages must be added to the specialty control on-chair package. Not available with Upholstered backs.

<input type="checkbox"/>	Quickie Sip n' Puff System		\$2,300.00
	116ID34	E2325 Sip N' Puff Kit	\$1,465.00
	116I01	E0955 Headrest	\$230.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
	116IH1	E2326 Sip n' Puff Tubing and Hardware	\$375.00
<input type="checkbox"/>	ASL Head Array Package (Adult)		\$4,140.00
	116I02	E2330 ASL Proximity Head Array (Adult)	\$3,910.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
<input type="checkbox"/>	ASL Head Array Package (Ped)		\$4,140.00
	116I03	E2330 ASL Proximity Head Array (Ped)	\$3,910.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
<input type="checkbox"/>	Switch-It Head Array Package (Adult)		\$4,140.00
	116I04	E2330 Switch-it Proximity Head Array	\$3,910.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
<input type="checkbox"/>	Switch-It Head Array Package (Pediatric)		\$4,140.00
	116I05	E2330 Switch-it Proximity Head Array	\$3,910.00
	116IJ01	E1028 Headrest Removable Mount	\$230.00
<input type="checkbox"/>	Quickie Proportional Chin Control Package		\$4,350.00
	116ID45	E2312 Proportional Joystick	\$4,120.00
	116IJ03	E1028 ASL Proportional S/A Joystick Mount	\$230.00

Attendant Control

<input type="checkbox"/>	116A02	E2331	Attendant Control for R-net	\$1,025.00
--------------------------	--------	-------	-----------------------------	------------

Attendant Mounting location

<input type="checkbox"/>	116DN1	-	Mounting location: Right Back Cane	N/C
<input type="checkbox"/>	116DN2	-	Mounting Location: Left Back Cane	N/C

Electronic Programmers

<input type="checkbox"/>	116PR5	-	R-net PC Programmer with Dongle	\$415.00
<input type="checkbox"/>	116PR6	-	DTT Rnet cables included.	\$1,030.00

11. POWERED SEATING

Powered Seating Control Method (TILT)

(For use with power seating actuators. Switch activation is stand-alone and does not require expandable or thru-drive electronics.)

<input type="checkbox"/>	116EO10	-	Switch Control 1 actuator w/dual toggle and mount. (For operating Tilt)	N/C
<input type="checkbox"/>	116EO11	-	Switch Control 1 actuator w/buddy button & mount. (For operating Tilt)	N/C

Switch Mount

(for use with switch control powered seating)

<input type="checkbox"/>	116EY01	-	Switch Mount Right	N/C
<input type="checkbox"/>	116EY02	-	Switch Mount Left	N/C

Powered Seating Thru-drive Control Method.

"Assignable Button" feature listed below (available on Rnet Systems) allows for direct "short-cut" access to various seating functions via joystick buttons. Functions accessed by depressing the button for more than 1 second. By selecting the thru-drive w/assignable button the chair will be pre-programmed with seating functions assigned to buttons on the joystick. Assignments can be changed via PC Dongle.

Thru-Drive 1 Actuator options.

<input type="checkbox"/>	116EO15	E2310	Thru Drive 1 Actuator w/out assignable buttons. For use with R-net (For operating Tilt thru Joystick.)	\$1,450.00
<input type="checkbox"/>	116EO13	E2310	Thru Drive 1 Actuator with assignable buttons. (For operating Tilt thru Joystick and/or Joystick. Requires R-net.)	\$1,450.00

12. WHEELS

Drive Wheel Type

<input type="checkbox"/>	116DW5	-	12" x 4" Pneumatic Tire	N/C
<input type="checkbox"/>	116DW6	-	12" x 4" Pneumatic w/Ins.	N/C
<input type="checkbox"/>	116DW27	-	12" x 3" Pneumatic Tire	N/C
<input type="checkbox"/>	116DW28	-	12" x 3" Pneumatic w/Ins.	N/C

The 12" x 3" wheel size reduces the overall width specification by approx. 2" and also alters the overall appearance of the chair.

Drive Wheel Tire

<input type="checkbox"/>	116DT1	-	Knobby Tire	N/C
<input type="checkbox"/>	116DT2	-	V Groove Tire (n/a with 12"x3" tires)	N/C

Wheel Lock

<input type="checkbox"/>	116WL26	-	Push to Lock	\$90.00
<input type="checkbox"/>	116WC1	-	Do Not Mount	N/C

Wheel Lock Extension Handle

Wheel-lock extensions handles are not available with 22" & 24" seats

<input type="checkbox"/>	116EH1	E0961	Extension Handle 6"	\$65.00
--------------------------	--------	-------	---------------------	---------

13. CASTERS

Front Caster Wheel

<input type="checkbox"/>	116CW17	-	8" Pneumatic	N/C
<input type="checkbox"/>	116CW20	-	8" Airless with Insert	N/C
<input type="checkbox"/>	116CW29	-	9" Pneumatic (not avail. With Suspension Fork)	\$90.00
<input type="checkbox"/>	116CW28	-	9" Pneumatic w/Insert (not avail. With Suspension Fork)	\$90.00

Caster Fork

Suspension Fork Options are only available with 8" casters.

<input type="checkbox"/>	116CF31	-	Solid Caster Fork	N/C
<input type="checkbox"/>	116CF35	E1016	Suspension Fork Med. 150-250 lbs. (Only avail. w/ 8" front casters)	\$270.00
<input type="checkbox"/>	116CF36	E1016	Suspension Fork Heavy. 251-350 lbs. (Only avail. w/ 8" front casters)	\$270.00

14. MANUAL CONVERSION KIT

<input type="checkbox"/>	116CV2	-	Manual Wheelchair Conversion Kit	\$1,025.00
--------------------------	--------	---	----------------------------------	------------

Axle

<input type="checkbox"/>	116AX1	-	Quick Release	N/C
<input type="checkbox"/>	116AX11	-	Quad Quick Release	\$95.00

Wheel Rim/Size

<input type="checkbox"/>	116RW1	-	24" Mag	N/C
<input type="checkbox"/>	116RW2	-	24" Spoke (not avail w/14" seat width)	N/C
<input type="checkbox"/>	116RW7	-	Omit Manual Rear Wheel	-\$55.00

Rear Wheel Tire

<input type="checkbox"/>	116RT1	-	Pneumatic	N/C
<input type="checkbox"/>	116RT2	-	Pneumatic Airless Insert	\$85.00
<input type="checkbox"/>	116RT3	-	Full Polyurethane	N/C
<input type="checkbox"/>	116RT4	-	Low Polyurethane	N/C
<input type="checkbox"/>	116RT6	-	High Pressure Clincher (not available with mag wheels)	\$100.00

Hand Rim

<input type="checkbox"/>	116HR1	-	Aluminum Anodized	N/C
<input type="checkbox"/>	116HR4	-	Plastic Coated	\$90.00
<input type="checkbox"/>	116HR7	-	Omit hand rim	-\$20.00

Projections

<input type="checkbox"/>	116HR8	-	8 Oblique	\$200.00
--------------------------	--------	---	-----------	----------

Wheel Locks

<input type="checkbox"/>	116WC2	-	Push to Lock Wheel Lock	\$45.00
<input type="checkbox"/>	116WC1	-	Do not mount	N/C
<input type="checkbox"/>	116ET1	E0961	6" Extension Handles	\$55.00

15. ACCESSORIES

Backpack

<input type="checkbox"/>	116A1	-	Black Utility Backpack	\$55.00
<input type="checkbox"/>	116A2	-	Black/Charcoal Kids	\$50.00

Seat Pouch

<input type="checkbox"/>	116A8	-	Black/Charcoal	\$50.00
--------------------------	-------	---	----------------	---------

Positioning Belt

<input type="checkbox"/>	116A13	-	2" Aircraft Buckle Padded	\$125.00
<input type="checkbox"/>	116A12	-	2" Aircraft Buckle	\$115.00
<input type="checkbox"/>	116A142	-	1.5" Auto Buckle	\$75.00

Wheel Chair Tray Table

Requires standard Armpads to mount.

<input type="checkbox"/>	116A31	E0950	Medium Tuffak (15"-17") Requires std. Armpads to mount tray	\$260.00
<input type="checkbox"/>	116A32	E0950	Large Tuffak (18"-20") Requires std. Armpads to mount tray	\$260.00
<input type="checkbox"/>	116A33	E0950	X-Large Tuffak (21"-22") Requires std. Armpads to mount tray	\$260.00

Miscellaneous

<input type="checkbox"/>	16A8105A1	-	Caddy	\$50.00
<input type="checkbox"/>	116A114	-	Dual Post Armrest Pouch	\$25.00

Additional Specifications

Weight:	115 lbs. (without batteries)
Overall Width	Using 12" x 4" Drive Tire: 24.5" using 14"-20" Seat Width 26.5" using 22" and 24" Seat Width
Overall Length:	31.5" (Without Footrests) Deduct 2" if using 12" x 3" Tire
Range:	As far as 30 miles (With Group 24 batteries depending on conditions and user environment)
Top Speed:	Up to 8.5 mph



Sunrise Medical (US) LLC
2842 Business Park Ave. · Fresno, CA 93727 · USA