



ORDER FORM

Q50 R Carbon

EIPW50

Customer	Ship To	
Account #	Name	
Date	Address	
PO#	City	
Buyer	State	Zip
Mark For	Phone	

Q50 R Carbon Totals	\$
Base Price	
Total of additional equipment	
Packaging and transport	
Total Price	

Base Specifications

User Weight Capacity	300 lbs.
Weight <i>(no batteries)</i>	32 lbs.
Seat Size	Width 17", Depth 17"
Drive Wheel Position	Rear-Wheel Drive
Base Size	Width 22.4", Length 36.3"
Drive Control	nVR2 Brushless Controller
Driving Range	7.5 miles (15 miles with additional battery)
Drive Motors	Brushless motor / 3.7 mph
Back Height	18"
STF Height	19.5"
Drive Wheels	8" with flat-free inserts
Tires / Casters	7" Solid
Charger	2A Off-Board
Battery Size	10Ah Lithium Ion
Color	Carbon Fiber with 5 color accents
Ground Clearance	2.75"

Additional Included Features and Accessories

Full Carbon Fiber

Quick-Release Folding

Lap Belt

Transport Strap

1.25" Thick Seat Cushion

Underseat Storage Bag

Storage for Additional Battery

Cup Holder

Off-Chair Battery Charging

Rear Running Lights

Depth Adjustable Joystick Mount

Flip-Up Footplate

Flip-Up Armrests

Breathable Back Upholstery

OPTIONS

OPTION	HCPC	DESCRIPTION	RETAIL
Model			\$
500MW93	K0899	Q50 R Carbon	\$3,595.00
Colors			\$
500C200	–	Black	NC
500C201	–	Red	NC
500C202	–	Blue	NC
500C203	–	Yellow	NC
500C204	–	Purple	NC
500C205	–	Orange	NC
Additional Battery			\$
500AG43	–	Additional 10 Ah Lithium Ion Battery	\$595.00

Date

Signature

Month

Day

Year

